

Zetasizer Nano Accessories Guide

MAN0487 Issue 2.0 June 2014

Copyright © 2007 - 2014 Malvern Instruments Ltd.

Malvern Instruments pursues a policy of continual improvement due to technical development. We therefore reserve the right to deviate from information, descriptions, and specifications in this publication without notice. Malvern Instruments shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this material. No reproduction or transmission of any part of this publication is allowed without the express written permission of Malvern Instruments Ltd.

Head office:

Malvern Instruments Ltd.
Enigma Business Park,
Groveswood Road,
Malvern,
Worcestershire WR14 1XZ
United Kingdom.

Tel + [44] (0)1684-892456

Fax + [44] (0)1684-892789

Zetasizer, Malvern and the 'hills' logo are registered trademarks in the UK and/or other countries, and are owned by Malvern Instruments Ltd.

NIBS and M3-PALS are trademarks of Malvern Instruments.
M3 is granted Euro Pat No: 1 154 266 DE FR.

Windows is a registered trademark of Microsoft Corporation.

Tygon is a registered trademark of Saint-Gobain Corporation.

Hellmanex is a registered trademark of Hellma GmbH & Co. KG.




Table of contents

Introduction and accessory range

Introduction	1-1
Accessory range.	1-1

General cells and cuvettes


Introduction	2-1
Cuvette holder	2-2
Cell and cuvettes	2-3
Size and molecular weight cuvettes	2-8
Folded capillary cell (DTS1070 / DTS1060/61)	2-12
Dip cell (ZEN1002)	2-20
High concentration cell (ZEN1010)	2-25
Surface zeta potential cell (ZEN1020)	2-29

Surface zeta potential cell

Introduction	3-1
Measurement technique	3-3
Preparation for measurement	3-4
Controlling an SZP measurement via an SOP	3-10
Performing the measurement.	3-16
Maintenance - cleaning the SZP cell	3-19
Chemical compatibility - SZP cell.	3-20

Flow-mode option

Introduction	4-1
Flow-mode	4-1
Applications.	4-2
Connectivity.	4-3
Exporting the flow-mode data	4-4
Editing and inspecting a flow-mode result	4-5



Controlling the flow-mode measurement via an SOP 4-9
The Flow-mode measurement display 4-14
Displaying the flow-mode report 4-15

Microrheology

Introduction 5-1
About DLS microrheology 5-2
Microrheology measurement process 5-3
Microrheology SOPs 5-4
Sample preparation advice 5-6
Running microrheology measurements 5-8
Analyzing microrheology data 5-13

Advanced protein features

Introduction 6-1
Controlling a protein mobility measurement via an SOP 6-4
Performing the measurement 6-7
Interpreting the results 6-8
Calculators tool 6-10

Introduction and accessory range

Introduction

This manual give an overview the accessories that are available for use with the Zetasizer Nano series of instruments.

This manual is a supplement to the following manuals:

- **Zetasizer Nano user manual**
- **Zetasizer Nano basic guide**



Warning!

The accessories or the samples to be measured may be hazardous if mis-used. Users must read the **health and safety** information in the **basic guide** before operating the system.

This manual focuses on specific issues of the Zetasizer Nano accessories that are not covered by the above manuals.

Accessory range

The accessories available for use with the Zetasizer Nano are indicted in the following section. Full descriptions on usage and application are contained in the subsequent chapters. The accessories that can be used will be dependent upon the instrument configuration and measurement type that will be performed.

Some accessories can be used for differing types of measurements, these are duplicated.

Cells and Cuvettes

Zeta potential measurements

DTS1070	Folded capillary cell Maintenance-free capillary cell primarily designed for zeta potential measurements. (This cell is a direct replacement for DTS1060/61).
ZEN1002	Dip cell Cell used to provide repeatable measurements of aqueous, and non-aqueous samples. It is particularly suitable for measurements of valuable aqueous samples where minimal sample quantity is important. Use with DTS0012 and PCS1115 cuvette - described below.
ZEN1010	High concentration cell Cell intended primarily for the measurement of zeta potential on a concentrated aqueous sample. The cell is suitable for a broad range of conductivities. It is particularly suitable for measurements of valuable aqueous samples where minimal sample quantity is important.
ZEN1020	Surface zeta potential cell Cell intended for the measurement of the zeta-potential at the surface of a flat material in an aqueous environment. Use with DTS0012 and PCS1115 cuvette - described below.

Size and molecular weight measurements

DTS1070	Folded capillary cell - This cell can be used for size measurements, in the Zetasizer Nano S, Nano ZS and Nano ZSP only. (This cell is a direct replacement for DTS1060/61).
DTS0012	Square polystyrene cuvettes - for size and molecular weight measurements.
ZEN0118	Disposable polystyrene low volume cuvette - for size and molecular weight (90° instruments only).
ZEN0040	Disposable plastic, micro cuvette, for size measurement at a 173° scattering angle.
PCS8501	Square glass cell with cap (round aperture) - for size and molecular weight. Also for use with the Dip cell.
PCS1115	Square glass cell with cap (square aperture) - for size and molecular weight. Also for use with the Dip cell.
ZEN2112	Low-volume quartz batch cuvette - for size and molecular weight.
ZEN0023	Quartz flow cell - for size, intensity measurements and molecular weight.

Note: All cuvettes have an outside dimension (o.d.) of 12mm

Cell to Zetasizer Nano instrument compatibility table

The table below indicates which cells and cuvettes are compatible with which instruments of the Zetasizer Nano range.

Cell	Zeta	Size	MW	S	Z	ZS	S90	ZS90	ZSP
DTS1070	•	•	•	•	•	•		•	•
ZEN1002	•				•	•		•	•
ZEN1010	•	•			•	•			•
ZEN1020	•	•			•	•		•	•
DTS0012		•		•		•	•	•	•
ZEN0118		•					•	•	
ZEN0040		•		•		•			•
PCS8501		•	•	•		•	•	•	•
PCS1115		•	•	•		•	•	•	•
ZEN2112		•	•	•		•	•	•	•
ZEN0023 (SEC systems)		•	•	•		•	•	•	•

The following chapters describe each cell or cuvette in more detail. The aspects covered are:

- Identification of each cell of the dispersion units with guidance on selection.
- How to fill each cell with sample and then insert into the Zetasizer Nano instrument.
- Some accessories require configuration of the software and SOP parameters. Where required the software controls are identified with explanation on how to use the cells to make measurements on the system.
- Maintenance procedures for inspecting and cleaning each cell and its respective components, including identification of the chemical compatibility of the cell components that may come into contact with the sample.

Additional accessories and options

A range of accessories and options are also available for more advance measurement strategies.

MPT-2 Titrator and Vacuum degasser

The MPT-2 Titrator is used to perform pH, additive and dilution titrations. The degasser is used to remove any dissolved gases from the titrants before dispensing by the Titrator. This prevents any gaseous escape during the titration increasing the accuracy and reproducibility.

Usage and operation is described in the **MPT-2 Titrator and Vacuum degasser user manual**.

- For use with **all** Zetasizer Nano instruments.

Flow-mode option

ZEN1006 Flow-mode option for Zetasizer Nano S and Zetasizer Nano ZS.

Instruments fitted with the flow-mode option can be connected to a size exclusion chromatography (**SEC**) system and be used as a light scattering detector.

The following kit is available for use with the Zetasizer Nano when connected to a SEC system. It includes the flowcell ZEN0023.

ZEN0116 Quartz flow cell kit for connection of the Zetasizer Nano to an SEC system.

Connection, use and operation of a flow-mode optioned Zetasizer Nano is described in the **Flow-mode** chapter later in this manual.

- For use with all Zetasizer Nano instruments **except** the Nano **Z**.
If fitted, an option part number label will be attached to the front of the cuvette holder. The flow mode option is included as a standard fitment on the Zetasizer Nano ZSP.

High Temperature

These instruments have an increased temperature range of 0° to 120°.

If fitted the Zetasizer Nano instrument label will include an '**HT**' identification and an option part number label will be attached to the front of the cuvette holder.

ZEN9063 Extends the upper temperature range of the Zetasizer Nano series from 90° to 120°

- For use with all Zetasizer Nano instruments **except** the Nano **ZSP**.

Narrow band filter

ZEN9051	Narrow band filter for 'Green' badged Zetasizer Nano S instruments
ZEN9052	Narrow band filter for 'Green' badged Zetasizer Nano ZS instruments
ZEN9061	Narrow band filter for 'Red' badged Zetasizer Nano S, Z and S90 instruments
ZEN9062	Narrow band filter for 'Red' badged Zetasizer Nano ZS, ZS90 and Nano ZSP instruments

This filter improves the signal for samples that fluoresce at the wavelength of the laser fitted. If a filter is fitted, an option part number label will be attached to the front of the cuvette holder.

- For use with **all** Zetasizer Nano instruments.

Microrheology

Microrheology is a new measurement type available to users of the Zetasizer Nano ZS and ZSP. It allows the measurement of the viscoelastic modulus of samples within the linear viscoelastic region.

Microrheology measurements require a software key to access the software features and functionality.

- For use with the Zetasizer Nano **ZS** and **ZSP**.

Protein mobility

Protein mobility is a new measurement type, supplied with the Zetasizer Nano ZSP as standard, and available to purchase separately for the Nano ZS.

The Zetasizer Nano ZSP has increased capability for the measurement of small and weakly scattering molecules and particles. This increased sensitivity improves the measurement of zeta potential, with its primary purpose being the ability to measure the zeta potential, or more appropriately the electrophoretic mobility of **protein** samples.

The Zetasizer Nano ZS is only capable of making these measurements at higher protein concentrations (>10-15 mg/ml).

A dedicated measurement type is included for protein mobility and also a suite of new calculators for proteins.

Protein mobility measurements are available after installing the **Advanced protein features** software key.

- For use with the Zetasizer Nano **ZSP and ZS**.

General cells and cuvettes

Introduction


Malvern offers a range of cells and cuvettes for performing measurements with the Zetasizer system. The choice of cell or cuvette is dependent upon the type of measurement being performed and the sample that will be measured.

The cells or cuvettes available for each measurement type are fully documented in this chapter with some discussion on their use and application. The aspects covered are:

- Identification of each cell of the dispersion units with guidance on selection.
- How to fill each cell with sample and then insert into the Zetasizer Nano instrument.
- Some accessories require configuration of the software and SOP parameters. Where required the software controls are identified with explanation on how to use the cells to make measurements on the system.
- Maintenance procedures for inspecting and cleaning each cell and its respective components, including identification of the chemical compatibility of the cell components that may come into contact with the sample.

Cuvette holder

The cuvette holder is for storing the cells before and after use. The cuvette holder swings out from under the instrument and up to 12 cuvettes can be stored.



ill 8503

**Caution!**

Ensure the thermal cap is lowered and all cuvettes have been removed before swinging the holder back under the instrument base.

The two trays that hold the cuvettes can be removed for cleaning.

The cuvette holder provides a place to store the **thermal cap** during changeover of cells. The cap is released by raising the cap and lifting off the cap post. Similarly, storage is also provided for the two “**thermal contact plates**” used with the folded capillary cell. When not being used, place these in the holder to the left of the tray.

The cuvette holder includes a **serial** number, **model** number and **option** labels. These identify the instrument and should be quoted in any correspondence with Malvern Instruments.

Cell and cuvettes

Which cell



Caution!

Due to the risk of melting, polystyrene cuvettes must not be used for measurements above 70°C.

The choices for each measurement type are outlined below with some discussion on their use.

Generally, for “easy to perform” measurements, such as with samples that scatter a reasonable amount of light (latex with 0.01% mass or higher, high scattering intensity, etc.) the disposable polystyrene cuvettes can be used.

- Disposable polystyrene cuvettes are easily scratched and should never be used more than once.
- Disposable cuvettes are not resistant to organic solvents, thus non-water based samples should generally be measured in glass or quartz type cuvettes.

The optical quality of the cells is vitally important when performing molecular weight and low concentration protein measurements (derived count rate <100kcps), therefore glass or quartz type cuvettes should be used to ensure the optimum signal is achieved.

Briefly the following cells can be used with the Zetasizer Nano instrument.

Cell	Application
Disposable “polystyrene” cuvettes – Standard and small volume	Size and zeta potential (with Dip cell)
Quartz glass cuvettes – Square, standard, low and ultra-low volume, flow	Size, molecular weight and zeta potential (with Dip cell)
Folded capillary cell	Size and zeta potential
High concentration cell	Size and zeta potential
Dip cell	Zeta potential
Surface zeta potential (SZP) cell	Zeta potential

All the cells mentioned are available from Malvern and should be used with the supplied cell caps. Using the caps will ensure greater thermal stability of the sample, as well as preventing dust introduction and possible spillage.

Cell and cuvette options

The cells and cuvettes described in this section can be used for all measurements.

	Folded capillary cell (DTS1070 / DTS1060/61)	Disposable polystyrene (DTS0012)
Application	Size, zeta potential	Size
Typical solvent	Water, water/alcohol	Water, water/ethanol
Optical quality	Good to very good	Good to very good
Minimum Sample volume	0.75ml	1ml
Advantages	Low cost Single use disposable (no cleaning) Use with MPT-2 Titrator No sample cross-contamination Fast sample change over (This cell is described later in this chapter)	Low cost Single use disposable (no cleaning) No sample cross-contamination
Disadvantages	Not resistant to organic solvents Unsuitable for use at high temperatures (above 70°C)	Not resistant to organic solvents Unsuitable for use at high temperatures (above 70°C)
Material	Polycarbonate	Polystyrene

	Disposable low volume polystyrene (ZEN0118)	Disposable low volume polystyrene (ZEN0040)
Application	Size	Size
Typical solvent	Water, water/alcohol	Water, water/alcohol
Optical quality	Good to very good	Good to very good
Minimum Sample volume	50µl	40µl
Advantages	Low cost Low volume Single use disposable (no cleaning)	Low cost Low volume Single use disposable (no cleaning)
Disadvantages	Requires careful filling to avoid bubbles Not resistant to organic solvents Unsuitable for use at high temperatures. (above 70°C) Only 90°C systems (Zetasizer Nano S90, ZS90, Zetasizer µV)	Requires careful filling to avoid bubbles Not resistant to organic solvents Unsuitable for use at high temperatures. (above 70°C) Only applicable to systems with NIBS optics (Zetasizer Nano S/ZS/ZSP)
Material	Polystyrene	Polystyrene

	Glass - square aperture (PCS1115)	Glass - round aperture (PCS8501)
Application	Size, molecular weight	Size, molecular weight
Typical solvent	Water, most organic and inorganic solvents	Water, most organic and inorganic solvents
Optical quality	Excellent	Excellent
Minimum Sample volume	1ml	1ml
Advantages	Highest optical quality Can use nearly any dispersant Reusable	Highest optical quality Can use nearly any dispersant
Disadvantages	Requires cleaning after measurement	Requires cleaning after measurement
Material	Glass	Glass

	Low volume Glass flow cell (ZEN0023)	Low volume quartz (ZEN2112)
Application	Size	Size
Typical solvent	Water, most organic and inorganic solvents	Water, most organic and inorganic solvents
Optical quality	Excellent	Excellent
Minimum Sample volume	75µl plus tubing	12µl
Advantages	Highest optical quality Can use nearly any solvent (tubing dependent) Use with MPT-2 Titrator	Highest optical quality Can use nearly any dispersant Low sample volume
Disadvantages	Requires cleaning after measurement With manual use requires careful filling to avoid bubbles	Requires cleaning after measurement Requires careful filling to avoid bubbles
Material	Glass	Quartz

	Low volume Glass flow cell pack (ZEN0116)
Application	Size as part of a Flow-mode (SEC) system Refer to ZEN0023 for description of cell Includes: - ZEN0023 flow cell - Connections for cell - 1.5m PEEK tubing (1/16 th inch outside diameter / 0.1mm inside diameter) Refer to the Flow-mode chapter later in this manual for connection, use and operation.

Use, cleanliness and filling advice

**Note**

Before filling and using a cell or cuvette, consult the **cleaning** section for each cell or cuvette, and perform any cleaning and maintenance procedures described.

When filling the cell there are several actions to consider; some that apply to all cells and others that are only applicable to the measurement type and the cell chosen.

The cleanliness of the cells used in each measurement is a paramount importance. As well as the information below, refer to the respective **cleaning cells/cuvettes section** for each cell, where more specific advice will be given.

Size cells and cuvettes

- All **size** cells should be rinsed/cleaned with filtered dispersant before use.

Molecular weight cells and cuvettes

- All **molecular weight** cells should be rinsed/cleaned with the filtered standard (e.g. Toluene) or solvent, then dried in a dust free environment such as a laminar flow cabinet, before use.

Zeta potential cells

- All **zeta potential** cells should be rinsed/cleaned with filtered dispersant before use.

Additionally

- The cell should be filled slowly to avoid air bubbles from being created. Ultrasonication can be used to remove air bubbles - but only if the sample is suitable for use with ultrasonics.
- If using syringe filters for the dispersant, never use the first few drops from the syringe, in case there are any residual dust particles in the filter that may contaminate the dispersant.

Size and molecular weight cuvettes


Filling a cell or cuvette

Fill the cell with the prepared sample as described below. Also refer to the filling advice given earlier in this chapter.

Standard cuvettes


A minimum sample volume must be provided. However, this minimum volume depends on the actual cell type and it is easier to ensure a certain **depth** of the sample in the cell.

This **minimum** is **10mm** from the bottom of the cell (the measurement is made 8mm from the bottom of the cell).



Do not overfill the cell, use about **15mm maximum**, as this can produce thermal gradients within the sample that will reduce the accuracy of the temperature control. ill 8506

- **When filling**, tilt the cuvette and allow it to fill slowly.
- To stop bubbles forming let the sample flow down the inside.

ill 7935

**Note**

When filled place a lid securely on the cuvette.

Low volume cuvettes

These cells are designed to use the minimum volume of sample possible for a size or molecular weight measurement. The sample must be pipetted carefully into the bottom of the cuvette, so it is filled from the bottom up.

The minimum volume that can be used is 12 microlitres for the Zetasizer Nano S, ZS and ZSP, 2 microlitres for the Zetasizer μ V and 20 microlitres for the Zetasizer Nano S90 and ZS90. This will only partly fill the visible cell volume. After filling, carefully inspect the cell for trapped bubbles.

**Note**

The lower sensitivity of the S90 and ZS90 means that they are unlikely to be suitable for the majority of molecular weight measurements.

Flowcells


Flowcells will be filled during the measurement procedure. Refer to the **MPT-2 Titrator and Vacuum degasser user manual**.

Inserting a cell or cuvette

For these measurements, perform the following:

Standard low volume cuvettes

- Open the cell area lid by pushing the button in front of the lid ①.
- Push the cell into the cell holder until it stops ②.
Some cells have opaque surfaces as well as polished optical surfaces. A polished optical surface must be facing the front of the instrument (towards the button). Most cells have a small **triangle** at the top to indicate the side that faces the front. This is especially critical for molecular weight measurements.
- Place the thermal cap over the cell ③ (this is not used if using a flowcell).
- Close the cell area lid ④.




ill 8510

Flowcells - using the MPT-2 Titrator

Follow these instructions for connecting a **flowcell** when using the MPT-2 Titrator. Always minimise the tubing within the cell area before inserting into the pinch valve channel.

- Follow the instructions in the previous section on how to insert the cell.
- Do **not** fit the thermal cap.



ill 8419

- The tubing is attached to the flowcell using threaded inserts; push the sample tubes into the inserts and screw into the top of the flowcell.
- The tubing is then inserted into the pinch valve channel; push **both** tubes down into the pinch valve on the side of the cell area. Ensure the join between the PTFE and silicone tubing is within the cell area.

Cleaning cuvettes

Two main types of cuvette are available:

Disposable polystyrene

Do not clean and re-use disposable cuvettes. It gives inaccurate results.


Reusable glass or quartz

The cleaning procedure depends on the sample measured so specific instructions cannot be given. Follow these guidelines:

- Rinse the cuvette with the dispersant that was used for the measurement.
- Try submerging the cuvette in an ultrasonic bath of clean solvent.
- Once clean, wipe the cuvette with a lint free tissue (photographers' lens cleaning tissues are recommended).
- The smaller and more dilute the sample, the more important cleanliness is.

Folded capillary cell (DTS1070 / DTS1060/61)

Description



(DTS1070 cell
shown above)

iii 8794

These are maintenance-free capillary cells primarily designed for zeta potential measurements, but **can** also be used for size measurements.

They have been designed to be used for a single measurement or series of measurements, then discarded rather than cleaned. This removes the chances of cross-contamination.

The cells are inserted with either the **Malvern logo** (DTS1070) or the **weld line** (DTS1060/61) facing the front of the instrument - refer to the **Inserting the cell** section later in this chapter.

The cells provide a low-cost alternative to previous reusable quartz capillary cells.

The stoppers can be replaced with 'Luer' connectors to provide leak-free connection to the optional MPT-2 Titrator.

Size measurements can also be performed without having to remove and reposition the cell.

Sample details can be written on the textured area on the side of the cell with a permanent pen.

Application	The cells are used for measurements of aqueous based samples
Typical solvent	Water, water/alcohol
Optical quality	Good to very good
Material	Body : Polycarbonate Electrodes : Gold plated beryllium/copper
Minimum Sample volume	0.75ml
Advantages	Low cost Single use disposable (no cleaning) Use with MPT-2 Titrator No sample cross-contamination Fast sample change over
Disadvantages	Not resistant to organic solvents Unsuitable for use at high temperatures (above 70°C)

Filling a folded capillary cell

The folded capillary cells should be filled with the prepared sample as described below. Both cells should be rinsed/cleaned with filtered dispersant before use; refer to **Cleaning the folded capillary cell** later in this section




Note

Filling a cell for a **protein mobility** measurement involves a different technique. Refer to the **Advanced protein features** chapter.

DTS1070 cell


The cell name - DTS1070 - will be identified on the central section of the cell body.

- Prepare the sample in a syringe of at least 1ml capacity.
- Place the sample syringe into one of the sample ports.
- Invert the cell ①.
- **Slowly** inject the sample from its syringe into the cell, filling the U tube to just over half way ②.
- Check no air bubbles form in the cell. Tap the cell gently to dislodge any that do form.



ill 7938

- Turn the cell upright and continue to inject slowly until the sample reaches the fill area as shown ③. Fill between **shoulder of cell** and the **FILL MAX** line.
- Check again for bubbles in the cell. Tap the cell gently to dislodge these.
- Check that the electrodes are completely immersed.



ill 8796 / 8797

- Remove the syringe and insert a cell stopper in **each** port.
- Remove any liquid spilt on the electrodes contacts.



Note


The stoppers **must** be fitted before a measurement is performed. Ensure that one stopper is fitted firmly, and the other one loosely, to avoid pressurisation of the cell.

DTS1060/61 cell

The cell name - DTS1060 - will be identified on the central section of the cell body.


- Prepare the sample in a syringe of at least 1ml capacity.
- Place the sample syringe into one of the sample ports.
- Invert the cell ①.
- **Slowly** inject the sample from its syringe into the cell, filling the U tube to just over half way ②.

- Check no air bubbles form in the cell. Tap the cell gently to dislodge any that do form.
- Turn the cell upright and continue to inject slowly until the sample is at the top of the electrodes ③.
- Check again for bubbles in the cell. Tap the cell gently to dislodge these.



- Check that the electrodes are completely immersed.
- Remove the syringe and insert a cell stopper in **each** port.
- Remove any liquid spilt on the electrode contacts.


iii 7938




**Note**

The stoppers **must** be fitted before a measurement is performed. Ensure that one stopper is fitted firmly, and the other one loosely, to avoid pressurisation of the cell.

Inserting the folded capillary cell



- Place a **thermal contact plate** into the recess on either side of the folded capillary cell ①. The plates provide increased temperature stability.



- Open the cell area lid by pushing the button ② in front of the lid.

1. Hold the cell near the top, away from the lower measurement area, and push into the cell holder until it stops ③.
2. The cells will fit in the cell either way round but do have a correct orientation. Please refer to the section below.
3. Close the cell area lid ④.


**Note**

Different versions of the capillary cells have unique thermal plates. If required, please contact your Malvern representative for the correct plates for your cell.

Folded capillary cell - orientation and insertion

The clear folded capillary zeta potential cells can show significant differences in sample count rate depending on the orientation of the cell in the cell holder. In most cases the difference does not affect the quality of the result, only the attenuator selection. However, in extreme cases, where the sample being measured is a poor scatterer, the measurement may not be possible in one of the orientations.

The diagram below shows the cell with the **preferred** orientation in the cell holder.



ill 8798 / 7945

DTS1070 cell

When inserting the cell, ensure that the **Malvern logo** faces towards the front of the instrument. Press down until the cell clicks into place.

DTS1060/61 cell

The cell is oriented such that the **weld line** is towards the front of the instrument. Press down until the cell clicks into place.


The cell is made of two different parts (front and back part), welded together.

Tests indicate that inserting the cell with the front part towards the laser gives better count rates, and hence this is the preferred cell orientation.

Using the MPT-2 Titrator

Follow the guidelines above on inserting the cell, then follow the instructions below for connecting the cell for the Titrator. Always minimise the tubing within the cell area before inserting into the pinch valve channel.

- The tubing is attached to the folded capillary cell using 'Luer lock' connectors.
- With a half-turn these secure to the Luer fittings on the top of the cell - do not overtighten.
- The tubing is then inserted into the pinch valve channel; push **both** tubes down into the pinch valve on the side of the cell area.



ill 8420

**Note**

The pinch valve manufacturer recommends that a vegetable-based oil (e.g. Castor oil) is used to lubricate the section of tube that is inserted into the pinch valve. This is done to help minimise friction, though testing by Malvern Instruments has not shown this to be essential.

Cleaning the folded capillary cell

This cell is intended to be used once then discarded. We recommend that, before a cell is used for the first time, it is flushed through with ethanol or methanol to facilitate wetting. A syringe or a wash bottle can be used. Use only sufficient fluid to wet the surface of the cell and electrodes.

The cell should then be flushed through with water as described below.


To clean the cell:

- Fill one syringe with de-ionised water or the dispersant.
- Place the full syringe in one of the sample ports on the cell and the empty syringe in the other.
- Flush the contents of the full syringe through the capillary cell into the empty syringe.
- Repeat the flushing process five more times, flushing the liquid backwards and forwards between the syringes. The cell is then ready for use.

Never attempt to clean the outside of the folded capillary cell. It causes small surface scratches that give inaccurate results.

Dip cell (ZEN1002)

Description



ill 8504

The Dip cell is used to provide a method to measure the zeta potential of both aqueous and non-aqueous samples. A number of samples can be prepared and the Dip cell inserted to measure each one in turn.

For aqueous samples the Dip cell can be used in conjunction with the disposable polystyrene (DTS0012). For non-aqueous samples use the reusable Glass - square aperture (PCS1115). These cells are described above.

Refer to **maintenance** and **chemical compatibility** section later in this chapter.

Application

The Dip cell can be used for measurements of aqueous and non-aqueous based samples.


Note

If the Dip cell is used for non-aqueous measurements, it is not recommended that it is subsequently used for aqueous measurements, as cleaning well enough afterwards to ensure the zeta potential standard is within specification cannot be guaranteed. It is recommended that two Dip cells are used: one Dip cell used for aqueous dispersants, and another for non-aqueous dispersants.

Filling the cuvettes used with the Dip cell

Fill the cell with the prepared sample as described below. Also refer to the filling advice given earlier in this chapter.

The Dip cell uses square cuvettes to hold the sample. With the insertion of the Dip cell the sample will be displaced upwards within the cuvette. If too much sample is placed into the cuvette prior to insertion of the Dip cell there is a risk that the cuvette will overflow.




ill 8505

To ensure a minimum sample volume is provided for the sample to be measured but protect against overfilling, we recommend the cuvette is filled to a depth of between **7mm** and **10mm** (**before** the Dip cell is inserted). The minimum level relates to approximately **0.7ml** of sample.

Do not overfill the cuvette; as well as overflowing the cuvette once the Dip cell is inserted, this can also produce thermal gradients within the sample that will reduce the accuracy of the temperature control.

- **When filling**, tilt the cuvette and allow it to fill slowly.
- To stop bubbles forming let the sample flow down the inside.




ill 7935

Once the Dip cell is inserted, It may be necessary tap the cell lightly to dislodge any bubbles that may be caught between the electrodes.

Inserting the Dip cell

Insertion of the Dip cell follows the same procedure as above, but first the Dip cell must be placed into the sample cuvette. This must be done at an angle to avoid any bubbles being caught between the sample electrodes.

- The cuvette **must not** be filled more than the recommended **maximum** depth of **10mm** ①.
- Tilt the cuvette to a maximum angle of 45° ②.
- Slowly insert the cell into the cuvette until the metal electrodes are covered ③. As the cell is inserted it displaces the sample so any bubbles will be pushed out from the top of the electrode gap.
- Once the electrodes are covered bring the cuvette up to the vertical ④.
- Inspect the combined cell and cuvette and check for any bubbles ⑤. If bubbles are present gently tap the bottom of the cuvette to dislodge these. If not dislodged repeat the above sequence.
- The cell can only be inserted one way round. Hold the base of the Dip cell cap and the top of the cuvette simultaneously ⑥. Ensure the coloured band on the label (and cuvette triangle) is facing the front of the instrument and push the cell into the cell holder until it stops - a '**stop**' on the Dip cell **must** rest on the top of the cell holder. Check that the cell sits flat on the cell holder.



ill 8508

**Note**


With the procedure complete, the measurement face of the cuvette (some have a small **triangle** at the top of the cell) and the coloured band on the Dip cell label must face in the same direction. This is to ensure the orientation is correct when inserted into the cell holder.

Removing the Dip cell

With care, by simultaneously holding the base of the Dip cell cap and the top of the cuvette, both the Dip cell and cuvette can be removed together. If adequate purchase cannot be obtained on both parts, then the following procedure is recommended.

- Lift the Dip cell up out of the cuvette, but before completely removing, gently tap the bottom of the Dip cell on the top of the cuvette ①. This will dislodge any remaining drops of sample from the cell into the cuvette.

If the Dip cell is simply lifted out of the cuvette there is a risk of drops of sample falling from the bottom of the Dip cell onto the instrument and surrounding area. This is especially important when using solvent based samples.



ill 8509

- Place the Dip cell immediately into an empty cuvette ② for storage. This will prevent any potential damage occurring either to the cell electrodes or the workspace.
- Remove the sample cuvette afterwards and place in the cuvette holder ③.

**Note**

Storage is also provided in the Dip cell case if the cell is not to be used for a while.

Cleaning the Dip cell


Clean the cell thoroughly between measurements, especially between different types of sample. Cross-contamination between samples can seriously affect the results.

**Caution!**

Do **not** immerse the complete cell. Only the sample electrodes must dip in to the dispersant, as shown below.

Clean the Dip cell electrodes regularly. They are made of solid palladium and can be cleaned physically and chemically. Follow the instructions below:

- Immerse the electrodes in a gentle ultrasound bath (30 Watts) for five to 15 minutes before use. Use the dispersant used for the previous sample as the cleaning fluid. If this dispersant contains additives such as surfactants, follow this by ultrasonicing for two minutes in the pure solvent.



ill 6763

**Warning!**

Take care: ultrasonication can produce a fine aerosol of the bath liquid.

- Remove the electrodes from the bath and rinse them with pure solvent. A pipe cleaner can be used for **gentle** cleaning of electrodes.
- To protect the Dip cell after cleaning, we recommend placing it in an empty cuvette for storage.

Before making a measurement, rinse the electrodes and cuvette with the sample to be measured.

When changing the sample, thoroughly rinse the electrodes with pure dispersant.

**Note**

The electrode holder is made from Natural PEEK (Polyetheretherketone) which is resistant to a wide range of chemical products. However, seek advice from Malvern and the sample manufacturer before using strong acid or base.


Chemical compatibility - Dip cell

With proper use, only the central electrode section of the Dip cell will ever come in contact with sample. The outer components of the Dip cell will only come into contact if spillage or overfilling occurs.

Component	Materials
Central section	
Electrode casing	Natural PEEK (Polyetheretherketone)
Electrodes	Palladium
Outer components	
Top and side casing	Natural PEEK (Polyetheretherketone)
Contacts	Phosphor Bronze with Nickel plating

High concentration cell (ZEN1010)

Description



ill 8450

The High concentration cell is intended primarily for the measurement of zeta potential of concentrated aqueous samples. The cell can be used in conjunction with the MPT-2 for automated titrations.

The cell consists of a high precision optical measurement block held within electrode chambers. This is all contained in an outer cuvette sized casing assembly that allows excellent thermal contact with the instrument cell holder.

Refer to **maintenance** and **chemical compatibility** section later in this chapter.

Application

The High concentration cell is used for measurements of high concentration aqueous samples.

The cell is supplied with the following components to prepare the sample and connect the cell:


1/32" internal bore silicon tubing with appropriate Luer fittings	Luer plugs for manual filling
Additional fittings for connection to the MPT-2 Titrator	Interdental brushes for cleaning of the electrode chamber, internal flow paths and optical block are also included

Filling the High concentration cell

Fill the cell with the prepared sample as described below. Also refer to the filling advice given earlier in this chapter.

Filling the High concentration cell uses a similar principle to the Folded capillary cell.


- Inject the sample **slowly** until the liquid reaches the bottom of the 'luer' outlet ①.
- Check no air bubbles form in the cell. Tap the cell gently to dislodge any that do form.
- Remove the syringe and insert a cell stopper in **each** port.
- Remove any liquid spilt on the electrodes.



ill 8448

Inserting the High concentration cell

The High concentration cell is inserted into the instrument and connected to the Titrator in the same manner as the Folded capillary cell.



ill 8449

The metal face of the cell must face the front of the instrument; this is to ensure good thermal contact between cell and instrument.

Cleaning the High concentration cell


General cleaning

Rinsing of the cell **prior** to a measurement should be carried out by flushing through with copious amounts of de-ionised water.

External surfaces of the assembled cell can be wiped clean with a weak soap solution.

Intensive cleaning

The cell first has to be disassembled before cleaning can be performed.



ill 8447

- Remove the screw cap ①.
- Separate the two halves of the cell ② by pulling the rear casing vertically away from the metal front.
- Note how the electrode chambers and quartz measurement cell block are assembled ③.
- Remove the chambers and cell block from the metal front casing ④.
- Detach the pipework and remove the top port ⑤.
- Protect the cell block from damage ⑥.

Once the cell has been disassembled, cleaning can be performed as described in the following table.

Component	Cleaning method
Screw cap	Wipe clean with a mild soap solution
Outer casing	Black part of casing (Rear - Delrin): Wipe clean with a mild soap solution. Metal part of casing (Front - Stainless steel): Immerse the casing in Hellmanex and place in a gentle ultrasound bath (30 Watts) for five to 15 minutes. Rinse with water once cleaned.
Electrode chambers and port	Electrode Chamber: Scrub gently with interdental brush and Hellmanex, then scrub with copious amounts of de-ionised water. Smaller internal bore: Scrub gently with interdental brush and Hellmanex, then scrub with copious amounts of de-ionised water.
Quartz measurement cell block	Scrub both internally and externally with interdental brush. Afterwards brush with copious amounts of water. Note: Once inserted back into assembly, a cotton bud with ethanol can be used for light cleaning of the outside of the cell block. This is only to remove any errant marks that may have occurred when assembling the cell.

Once cleaned, leave all parts to dry before re-assembling. Re-assembly is the reverse of dis-assembly. Take care not to damage the sprung electrodes located in the rear casing.


Chemical compatibility

With proper use, only the central electrode and measurement section of the High concentration cell will ever come in contact with sample. The outer components of the cell will only come into contact if spillage or overfilling occurs.

Component	Materials
Central section	
Electrode chambers / O-rings	Natural PEEK / Nitrile rubber
Electrodes	Palladium
Electrode contacts	Brass
Precision measurement block	Quartz
Tubing	Silicone rubber
Outer components	
Casing	Delrin / Stainless steel 316
Cap	Delrin / Phosphor Bronze
Contacts	Gold plated beryllium / Copper

Surface zeta potential cell (ZEN1020)

Description



ill 8687

The surface zeta potential cell is intended for the measurement of the zeta potential at the surface of a flat material in an aqueous environment. The cell is a dip cell type device and can be used with 1ml cell DTS0012 and PCS1115. It is **incompatible** with the MPT-2 Titrator.

Application

The surface zeta potential cell is used for measurements of aqueous samples.


For details for this cell, refer to the separate **Surface zeta potential** chapter in this manual.

Surface zeta potential cell


Introduction

This chapter gives an overview of the Zetasizer Nano cell for measuring surface zeta potential. It describes how to use, insert and clean the cell to ensure reliable and consistent measurements.

The **Surface zeta potential (SZP)** cell is intended for the measurement of the zeta potential at the surface of a flat material in an aqueous environment. The cell is a dip cell type device and can be used with 1ml cell DTS0012 and PCS1115. It is **incompatible** with the MPT-2 Titrator.



The cell consists of a height adjustable sample barrel ⑤, in which the sample is glued onto a sample holder ⑥ and held between two palladium electrodes ⑦. A series of zeta potential measurements are then performed in a conventional cuvette ③, with the measurement position within the cell controlled by rotating the cell cap ① which adjusts the height of the sample barrel.



The cell is supplied with the following components to prepare, load and set the sample:


Surface zeta potential cell with palladium electrodes	A 12-well plate for storing the samples
10 PEEK sample holders	A screwdriver for cell tightening
Forceps for sample handling	A cell height alignment tool and a sample holder for gluing the sample to

Measurement technique

A surface zeta potential measurement consists of attaching a sample to a mount or holder that is then held in place between two electrodes. The sample is then immersed in an appropriate aqueous solution, containing tracer particles.

The apparent tracer mobility is now measured at a number of different distances from the sample surface. The electro-osmotic flow at the sample surface will tend to fall off with increasing distance hence; close to the surface the tracer mobility will be dominated by the electro-osmotic surface flow, while at distances further from the surface it will be dominated by the electrophoretic motion of the tracer itself.

The graph below shows a typical plot of reported zeta displacement from the surface.



The zeta potential at the surface is then calculated by extrapolating the graph to zero displacement and applying the following formula; ill 8694

$$\text{Surface zeta potential} = -\text{intercept} + \text{tracer zeta potential}$$

where the tracer zeta potential is recorded far from the wall, where the electro-osmotic flow can be taken as zero.

In the displacement graph above, the **blue circles** represent the reported zeta potential of the tracer particles, while the **red squares** represent the zeta potential of the tracer particles measured far from the sample surface and also independent from any electro-osmotic effects.

Preparation for measurement

Before a measurement can be performed the cell must first be loaded onto the sample holder and then attached to the cell. The complete cell is then inserted into a standard cuvette and placed into the instrument.

These following operations are described in the next sections; once these are complete the measurement can be performed:

- Loading the surface zeta potential cell with sample.
- Inserting the surface zeta potential cell into the instrument.


Loading the surface zeta potential cell

The surface zeta potential cell is loaded with a sample as described below:

**Note**

Take care not to damage the sample surface during attachment to the sample holder.

The sample to be measured should be cut into rectangular pieces no larger than 7mm x 4mm (LxW) and no more than 1.5mm thick (H).



ill 8689

The sample is then attached to the sample holder using an appropriate adhesive such as Araldite™ (refer to note overleaf for guidelines on glue selection).

- A sample gluing tool is provided to hold the sample holder during the gluing process.

Ensure the sample is placed squarely onto the sample holder and not at an angle. The sample should be perpendicular to the electrodes, once inserted into the barrel. There should be no large gaps between the electrodes and sample ($\leq 200\mu\text{m}$).

Wait until the sample has set in place, then load onto the cell using the supplied forceps and screw the holder into place.

**Note**

The glue used should be selected beforehand and be compatible with the experimental design. It should be capable of attaching to both the sample and the holder, and it should not be soluble in the selected medium, so that the sample is securely held in place for the duration of the experiment.

Inserting the surface zeta potential cell

The insertion of the surface zeta potential cell into the Zetasizer Nano is done in three stages.


- First a **coarse** alignment is performed where the cell has to be aligned to a zero position with respect to the instrument laser - this is initially done using a height alignment tool supplied with the cell.
- Secondly the cell is inserted into the cuvette, which is then added into the instrument.
- Thirdly a **fine** alignment is performed using the count rate meter in the application software.

Once these stages are complete a measurement can be performed.


Zeroing the cell height position - coarse alignment

Once the sample holder is in place, the surface of the sample must be aligned to a zero height position with respect to the instrument laser, using the height alignment tool.

The surface is aligned to a **zeroing** target etched on the windows of the tool. There are two **zeroing** targets, one on the front plate, and one on the back. Hold the alignment tool so that the centre of the two targets coincide, to avoid a parallax error.



Insert the cell assembly into the tool, so that the white mark on the cell is facing the front of the tool, indicated by the **white spot**, and tilting forward. Adjust the cell cap to alter the sample barrel position until the surface of the sample is aligned with the zeroing target on the tool window. ill 8690




ill 8691

The sample barrel position should be adjusted so that the sample surface and the centres of the two zeroing targets all line-up exactly.

With the sample height set, the cell cap also needs to be zeroed.

Loosen the cap screw, then rotate the cell cap until the white mark on the cap is in line with the white mark on the cell body. Secure the screw afterwards.



ill 8692

Filling the cuvette and inserting the cell into the instrument

The surface zeta potential cell can use either disposable plastic and glass or quartz square cuvettes to hold the sample. With this insertion of the cell, the dispersant will be displaced upwards within the cuvette.

- Prepare an appropriate aqueous suspension, containing the tracer particles. Sufficient suspension should be added to the cell so that the sample, electrodes and screw are all completely submerged; This is approximately 1.2 ml.
- Fill the cuvette the prepared aqueous suspension. When filling, tilt the cuvette and allow it to fill slowly.
- To stop bubbles forming, let the sample flow down the inside.


With cell filled, the cell can be inserted; this must be done at an angle to avoid any bubbles being caught between the sample electrodes, and to ensure that the sample plate is entirely submerged. The level of the dispersant must be significantly above the top of the electrodes and the nylon screw that holds the sample in place. The procedure follows below.



Note

With the procedure complete, the measurement face of the cuvette (some have a small triangle at the top of the cell) and the white mark on the cell body must face in the same direction. This is to ensure the orientation is correct when inserted into the cell holder.

- The cuvette must not be filled more than the recommended maximum depth of 20mm before insertion of the cell ①.
- Tilt the cuvette to a maximum angle of 45° ②. This is to avoid spilling the dispersant.
- Slowly insert the cell into the cuvette until the sample holder, barrel and electrodes are covered ③. As the cell is inserted it displaces the sample so any bubbles will be pushed out from the top of the electrode gap.
- Once the electrodes are covered bring the cuvette up to the vertical ④.
- Inspect the combined cell and cuvette and check for any bubbles ⑤. If bubbles are present around the electrode or nylon screw, gently tap the bottom of the cuvette to dislodge these. If not dislodged repeat the above sequence.



III 8799

- Hold the base of the SZP cell cap and the top of the cuvette simultaneously ⑥. Ensure the white marks ① on the cap and cell body, and the cuvette triangle ②, are facing the front of the instrument and push the cell into the cell holder until it stops - a 'stop' on the surface potential cell must rest on the top of the cell holder. Check that the cell is sitting flat, and that the cuvette is also **fully** inserted and rests on the base of the cell holder.

Fine alignment using the Count rate meter

A final fine adjustment of the zero-position can be made once the cell has been placed into the instrument; this is done using the **Count rate meter** in the Zetasizer software.

From the main menu select **Tools-Count rate meter** to open the count rate meter, and set the count rate meter as specified:

- Select the **Forward scatter** radio button.
- Under **Cell type**, select **ZEN1020 plate cell** (Surface zeta potential cell) from the drop down list.
- Set the **Attenuator** to **11**.

To identify if your sample is aligned correctly open the instrument lid and rotate the cell cap **clockwise** in increments of approximately 1/8 to 1/4 of a turn. Close the lid and observe the count rate measured in between each increment. Once you lower the sample far enough that the **count rate** observed is **zero**, open the lid and rotate the cell cap **counter-clockwise** by approximately 1/8 to 1/4 of a turn. Close the lid and observe that the **count rate** has **risen** a reasonable extent. This position is the starting point for your measurement.

For ease of reference, it is now possible to loosen the thumb screw on the micrometer and rotate it to the front of the cell then retighten it to proceed with measurement.



Note

After the fine alignment step, whenever the cell height is adjusted during the experiment, it is imperative that the physical position of the cell is not moved within the cell holder. Any movement of the cell will result in a different zero height and this must be constant throughout a given measurement. Any alteration of the height during an experiment will reduce the quality of the data.

Controlling an SZP measurement via an SOP

Once the cell has been inserted into the Zetasizer and zeroed, a measurement can be made. This is done using an SOP, or a manual measurement in the usual manner.

An SOP can be configured to control all settings for the accessories automatically.

A **surface zeta potential (SZP)** measurement follows the same SOP format as performed when doing a normal zeta potential measurement, with a few exceptions. When a SZP measurement is chosen two extra dialogues - **SZP measurement** and **Tracer measurement** - will be included in the SOP selections.

The **SOP Editor** and setup is described in full in the **Zetasizer Nano user manual**. Most of the SOP sections are common to Measurement types, and these are described in the above manual. The other SOP sections are specific to the SZP cell being used; these are described below. Also note that some of the other dialogue pages will alter slightly to accommodate extra parameters necessary to perform the SZP measurement.

Creating or editing an SOP - Measurement Type selection

- To create a new SOP, select **File-New SOP**. This will open up the SOP Editor. The SOP Editor consists of several dialogues that can be stepped through by using the **Next** arrow button.
(To edit an existing SOP, choose **Open-SOP** instead.)
- Complete the SOP Editor as described in the Zetasizer Nano user manual.
- Once the SOP has been created, press **Finish** and save the new SOP.

The various SOP dialogues are described below.

Measurement type options


Select a **surface zeta potential** measurement type then complete the SOP creation as required.

Tracer Material

Please refer to the **Sample - Material** description in the **zeta potential SOPs** section of the main user manual.

Cell

Please refer to the **Sample - Cell** description in the **size SOPs** section of the main user manual.




As the surface zeta potential cell was selected as the measurement type, this cell will be the only cell choice available. The Zetasizer software will configure all settings and parameters to match this cell.

The default selection is the **Surface zeta potential cell**.

SZP measurement

The **SZP measurement** SOP window is similar to the standard zeta potential Measurement window. Where appropriate please refer to the **measurement** description in the **zeta potential SOPs** section of the main user manual for more details on each of the measurement options.




Settings	Description
SZP measurement duration	The SZP (Surface zeta potential) measurement duration options are the same as standard measurement duration options available during normal zeta potential measurements.
SZP measurements	The SZP measurements options define the number of repeat measurements made at each displacement away from the surface and the length of any delay between repeat measurements.
SZP displacement	The SZP displacement options define the distances to be used during a surface zeta potential measurement. The Number of positions defines the number of points away from the surface where tracer mobility is to be measured. The Size of steps defines the additional distance away from the surface that each measurement is made. Note that 500 microns can only be chosen if the tracer measurement displacement (in the tracer measurement SOP window) is set to 1500 microns or greater. This is to ensure that a minimum of three displacement points will be measured.

SZP measurement - Advanced

Refer to the **Measurement - Material** description in the **zeta potential SOPs** section of the main user manual.

Tracer measurement

The **Tracer measurement** SOP window is similar to the standard zeta potential Measurement window. Where appropriate please refer to the **Measurement** description in the **zeta potential SOPs** section of the main user manual for more details on each of the measurement options.



Settings	Description
Tracer measurement duration	The Tracer measurement duration options are the same as standard measurement duration options available during normal zeta potential measurements.
Tracer measurements	The tracer measurements options define the number of repeat measurements that can be made at each measurement distance from the sample surface, and the length of any delay between these repeat measurements.
SZP displacement	The final stage in a surface zeta potential measurement is a fast field reversal (FFR) only measurement. The purpose of this is to make a measurement only of the tracer mobility, which will not include any electro-osmotic component, and this will be used in the surface zeta potential equation.


Settings	Description
SZP displacement (continued)	The Tracer measurement displacement defines the distance from the sample at which this FFR only measurement takes place. The displacement is altered in 125micron increments.

Tracer measurement - Advanced

Refer to the **Measurement - Material** description in the **zeta potential SOPs** section of the main user manual.

Data processing

This window allows the advanced analysis parameters to be set. It is generally best to leave these set to default.



Size ranges and measurement thresholds can be applied to the analysis to filter spurious peaks prior to the analysis being performed. These can be setup using the **Configure** button.


Settings	Description
Configure button	<p>Pressing the Configure button will display the Plate analysis parameters window, which enables various attributes of the analysis model to be altered. These include the measured zeta Display range, and the measurement thresholds.</p> <p>If it is known that all particles within the sample will fall within a certain zeta potential range, then the zeta Display range can be set to improve the repeatability of the measurement result; similarly a lower threshold sets the noise rejection baseline in the zeta potential distribution.</p>

Refer to the **Help** file for more information.


Performing the measurement

With the cell loaded into the instrument, and the SOP configured, a measurement can be performed.

When the measurement is started, a user instruction is given to turn the cap on the top of the cell by a given amount - this will set the distance to the first required displacement position. This is a manual operation and the user must open the cell area lid, turn the head of the cell the specified amount, then close the lid again before continuing.



Each 1/4 turn of the cap counter-clockwise corresponds to a movement or displacement of 125 microns; a second 1/4 turn will correspond to a total displacement of 250 microns, and so on. The amount to move the cap will be indicated in the SOP and in the measurement instructions.



Once all of the specified measurements have been performed at that displacement, an instruction is given to set the cell to the next measurement position. This process will continue until all measurements have been made at all the positions specified in the SOP, and the surface zeta potential measurement is then complete.

The data is stored as a 'parent' surface zeta potential record, with 'child' records relating to the individual zeta potential measurements made at each displacement.

Record	Type	Sample Name	Measurement Date and Time	T °C	ZP mV	Mob µcm/Vs	Cond mS/cm
139	Surface Zeta		28 June 2011 14:17:10	25.0			
140	Surface Zeta		28 June 2011 14:17:36	25.0			
141	Zeta	1 125 microns 1	28 June 2011 14:17:42	25.0	-76.2	-5.977	0.279
142	Zeta	1 125 microns 2	28 June 2011 14:18:03	25.0	-76.2	-5.977	0.279
143	Zeta	1 125 microns 3	28 June 2011 14:18:09	25.0	-76.2	-5.977	0.279
144	Zeta	2 250 microns 4	28 June 2011 14:18:24	25.0	-72.5	-5.683	0.279
145	Zeta	2 250 microns 5	28 June 2011 14:18:31	25.0	-72.5	-5.683	0.279
146	Zeta	2 250 microns 6	28 June 2011 14:18:37	25.0	-72.5	-5.683	0.279
147	Zeta	3 375 microns 7	28 June 2011 14:19:02	25.0	-68.7	-5.389	0.279
148	Zeta	3 375 microns 8	28 June 2011 14:19:08	25.0	-68.7	-5.389	0.279
149	Zeta	3 375 microns 9	28 June 2011 14:19:15	25.0	-68.7	-5.389	0.279
150	Zeta	4 500 microns 10	28 June 2011 14:19:24	25.0	-65.0	-5.095	0.279
151	Zeta	4 500 microns 11	28 June 2011 14:19:30	25.0	-65.0	-5.095	0.279
152	Zeta	4 500 microns 12	28 June 2011 14:19:37	25.0	-65.0	-5.095	0.279
153	Zeta	0 0 microns 13	28 June 2011 14:19:49	25.0	-50.1	-3.931	0.279
154	Zeta	0 0 microns 14	28 June 2011 14:20:08	25.0	-50.1	-3.931	0.279
155	Zeta	0 0 microns 15	28 June 2011 14:20:13	25.0	-50.1	-3.931	0.279
156	Surface Zeta		28 June 2011 15:42:59	25.0			
Mean 140				25.0			
Std Dev							

A surface zeta potential report is available to view the results. Select **View-Workspaces-Surface zeta potential** to view the appropriate workspace.

Editing the results

Surface zeta potential results can be edited by right-clicking on the record in the **records** view and selecting **Edit result**; or select **Edit-Edit result** from the main menu. With the **Edit result** window open, the Debye length model, tracer material and dispersant properties can then all be changed.

The surface zeta potential edit result option allows points to be removed from the displacement plot by left-clicking on them on the displacement graph. When **OK** is clicked, a new surface zeta potential record is created containing only the child zeta potential measurements that were included in the analysis.

Removing the cell from the instrument

Follow the procedure for removing the **Dip cell** from the instrument; the operation is the same.

Maintenance - cleaning the SZP cell

**Caution!**

During cleaning and use it is vital not to let any fluid enter the top and cap area of the cell assembly.

Any cross contamination of material from one measurement to the next could affect the result, so it is extremely important to ensure the cell is completely clean before use.

Cuvettes used with the surface zeta potential cell


- If a quartz cuvette was used for the measurement, it is recommended to clean the cuvette with Hellmanex, and then rinse with copious amounts of de-ionised water, **prior** to reusing it.
- If a plastic disposable cell was used for the measurement, it is recommended that this is disposed of and a new one used for all subsequent measurements.

General cleaning

As a complete assembly the cell can be cleaned using de-ionised water or with a Hellmanex solution. If Hellmanex is used, the cell **must** be rinsed with copious amounts of de-ionised water, **prior** to reusing it.

More efficient cleaning can be obtained by immersing the electrode area and sample holder in a gentle ultrasound bath (30 Watts) for 5 to 15 minutes.

- Use the dispersant used for the previous sample as the cleaning fluid. If this dispersant contains additives such as surfactants, follow this by ultrasonicing for two minutes in the pure solvent.



Over time, it is likely that the electrodes will become discoloured or tarnished. This is expected, and although it cannot be cleaned, this will not affect the quality of the data obtained. ill 6763

Intensive cleaning

Cleaning can be performed as described in the following table. The material and chemical compatibility of each component is detailed in the next section.

Component	Cleaning method
Cell cap	Wipe clean with a mild soap solution. Rinse with water once cleaned.
Outer casing Sample barrel Sample holder	Wipe clean with a mild soap solution. Rinse with water once cleaned.
Electrodes	Scrub gently with a pipecleaner and Hellmanex, then scrub with copious amounts of de-ionised water.

Once cleaned, leave all parts to be fully dry before re-using; especially the electrode and sample holder area.

Chemical compatibility - SZP cell

Components of the Zetasizer Nano that may come into contact with the sample are manufactured from materials that are considered to give the widest protection from chemical attack. However, it is important to check that any sample or titrant used is chemically compatible with the materials mentioned.



Warning!

It is advisable that the chemical compatibility is checked against the materials identified below before inserting a sample. It is also recommended that a test is performed on the material with the sample before more permanent usage is undertaken.

With proper use, only the central measurement section (see table for components) of the surface zeta potential cell will ever come into contact with sample. The outer components of the cell will only come into contact if spillage or overfilling occurs.

Component	Materials
Casing	Natural PEEK / Stainless steel 316
Sample barrel	Natural PEEK
Sample holder	Natural PEEK
Sample holder screw	Nylon
Electrodes	Palladium
Contacts	Beryllium / copper

Flow-mode option

Introduction

This chapter gives an overview of the **Flow-mode** option. Instruments fitted with this option can be connected to a size exclusion chromatography (**SEC**) system and be used as a light scattering detector.

ZEN1006 Flow-mode option for Zetasizer Nano S and Zetasizer Nano ZS.

It describes how to connect, control and operate the flow-mode arrangement to ensure reliable and consistent measurements.

The following cell is available for use with the Zetasizer Nano when connected to a SEC system.


ZEN0116 Quartz flow cell kit for connection of the Zetasizer Nano to an SEC system.

- The option is for use with all Zetasizer Nano instruments **except** the Nano **Z**. If fitted, an option part number label will be attached to the front of the cuvette holder. The flow mode option is included as a standard fitment on the Zetasizer Nano ZSP.

Flow-mode

The **Flow-mode** option allows the Zetasizer Nano to be connected, using a flowcell, to a flowing sample stream, such as the output of a chromatographic column or a field flow fractionation system, and measurements made without interrupting the flow. The output from external detectors, such as refractive index and ultra-violet absorption detectors, can be input back into the Zetasizer to allow integration of the data, using optional hardware.

In a flow-mode measurement, both the scattered light intensity and hydrodynamic diameter are plotted as a trend, and the addition of optional analogue inputs enables the simultaneous display of data from up to two other detector outputs.



ill 7908

By connecting the output from the external device to the external input socket on the **rear** of the Zetasizer Nano, a real-time parameter reading from the device can be directly inputted back into the Zetasizer Nano software. This parameter reading can be plotted as a trend, thus enabling additional sample characteristics to be monitored.

Applications

Applications for this feature include use as a chromatography detector and a process monitor.


Separations detector

The Zetasizer Nano can be connected to a size-exclusion chromatography (SEC) system or a field-flow fractionation (FFF) system.

When connected to one of these systems, the Zetasizer Nano acts as a sensitive light scattering detector, simultaneously plotting trends of the total intensity of light scattered and the size as the hydrodynamic diameter.

Light scattering is an almost universal detector, as most materials, such as proteins for example, will have a different refractive index to the buffer they are in, and hence will produce a light scattering signal.

Other detectors can be connected to the Zetasizer Nano optics unit and their signals plotted on the same axes as the light scattering signal. The timing of these signals can then be adjusted to compensate for the output delay due to the detectors being connected serially in the flow path. These detectors will usually be refractive index or UV and used to measure the concentration of the eluting sample. These values can then be used to calculate the molecular weight of each component separated.



ill 7906

Process monitor

The Zetasizer Nano can be connected to a flowing sample extracted from a process stream or reaction vessel. The high concentration capability of the Nano S and ZS means that many processes can be monitored without further sample preparation, simplifying the measurement.




Note

The method of sampling the process and transferring this to the Zetasizer is not provided.

Connectivity

To obtain the data from the external detector use the 4-way 'Lemo' socket on the rear panel.

A lead is supplied - the end terminated with a 'Lemo' plug connects to the rear panel; whilst the other end, terminated in 4 bare wires, will be connected to the external device. External devices can be connected to the external input 1 and earth, and external input 2 and earth. An additional input may be used as a trigger for remotely starting the Zetasizer Nano.



The lead supplied by Malvern can be connected directly to the rear of the external device, or to an output from the device, whichever is appropriate.

Consult the external device documentation to obtain the output connections for connecting the Malvern lead.

Connector input voltage specification

The input voltage is -5v to +5v (analogue).

Exporting the flow-mode data

Once received, the external device data obtained during the flow-mode measurement can be exported from the Zetasizer Nano software, saved as a text file, and inserted into a spreadsheet software package (such as Microsoft Excel) for analysis.

This is done using the **Export flow result** macro option in the **Tools** menu.


If this macro is not visible in the menu, select **Tools-Options...** to enable it.

Exporting the flow-mode data:

- Select a flow measurement record. The export will not work if this is not done.
- Select **Tools-Macros-Export flow result**.
- A window is displayed requesting the export destination.
- On selecting **Save**, the data will be exported and saved as a text (.txt) file. The file can then be imported into the target program.

Editing and inspecting a flow-mode result

A flow-mode result can be edited by right-clicking on the measurement record and selecting **Edit-Result**. This will display the window below.





ill 8822

The window shows the following:

① Chromatogram

This shows the flow-mode plot. By default the **Intensity** trace is displayed on the left y-axis, and the **Z-average mean** parameter displayed on the right y-axis. Right-clicking on the y-axis text allows other available traces to be displayed. The available traces are displayed in the **Selected trace** area of the **Traces** tab.

-
- Ⓐ **Traces**
- The traces on the plot display the measured parameters with respect to the volume passed through the flow-mode system arrangement. Use the cursor to select each trace; the trace selected will be displayed in the plot legend and the displayed **Trace** tab.
- The traces shown can be chosen by right-clicking on the y-axes or by using the **Trace** tab.
- Generally the **Z-average mean** values are shown as dots on the plot, whilst the **Intensity** is shown as a continuous line.
-
- Ⓑ **Peaks**
- The **peaks** detected by the analysis are indicated by the **blue vertical bars** on the plot. The peaks displayed are determined by inspecting the Z-average mean data points.
- When a peak is selected, by either clicking on the bar or selecting the peak check box, the **Peak** tab will be displayed. This shows the peak details (see below). Note that the selected peak will be highlighted deep blue.
- By clicking and dragging, it is possible to select a region of the chromatogram, which is highlighted in purple. By right clicking the mouse, it is then possible to perform a DLS analysis of this region. This calculates and intensity-weighted average Z-average that is calculated in the box in the bottom right of the display.
- Note that when the user clicks **OK**, although the data is saved, it is not displayed in the reports and must be accessed by returning to the **Edit results** dialogue.
-
- Ⓒ **Viewing tools**
- These enable different plot views to be shown. The tools can be selected and used in two ways:
- Select directly from the drop-down menu and perform the action.
 - With the drop-down menu showing the Peak/trace highlighter, hold the appropriate control key down to perform the action.
- The actions available are:
- Peak / Trace highlighter (cursor)** 
- When the cursor is moved over the traces, peaks and axes of the plot, this will change to a hand, allowing the feature to be selected. When a trace or peak is selected the respective tab will be displayed.
- Zoom**  (Alt)
- Use to place a “marquee” around a specific part of the plot and enlarge. The axes will automatically adjust to match the enlarged plot.

 © Viewing tools
(continued)

Scale  (Control)

Use to dynamically zoom in or out of the plot. This will be centred on the cursor position. The axes will automatically adjust to match the new scale.

Pan  (Shift)

Allows the plot to be moved left, right, up and down. The axes will automatically adjust to match the new position.

 Reset view

Resets the view back to the initial view. All zoom, scale and pan actions will be removed.

 Ⓓ Baseline
pointer

The **Baseline pointer** is used to set the datum level for the plot. When the measurement is originally performed the software will attribute an appropriate baseline level and use this as a filter to remove unwanted and erroneous data. All measurement data below this level will be removed from the measurement, plus any data above it that the software considers to be uncharacteristic of the measurement, in effect an “outlier”.

This filter of outliers can be disabled in the **Data Processing** section of the SOP and **Edit results** dialogue. Choose **Configure** and then disable the option to **Filter flow mode data** to do this.

With reference to the originally attributed baseline level, the software will display any peaks that may have occurred during the measurement.

To observe the effect on the peaks when the baseline is adjusted, change a y-axis to **Intensity**, select the baseline pointer with the cursor and drag to a different position along the right y-axis. The peaks will be displayed in the **Peaks** tab (see below) - this will automatically be displayed when the baseline pointer is moved.

Note that when the baseline level is moved the software will remove any data that is below the new level. This may significantly alter any peaks shown.

The baseline pointer is only visible when an **Intensity** trace is displayed.

 Ⓔ Trace data

Sections of the chromatogram can be removed from the flow-mode plot. To remove a part of the chromatogram, highlight a section by dragging the mouse over the graph; then right-click on the highlighted section and select **Trace data**.

Options will be given to either remove the highlighted section (**Delete data in this section only**), or to remove all data outside the highlighted section (**Keep data in this section only**).

② Measurement tab

This **Measurement** tab displays the measurement details. These are the same details as entered when the measurement was first performed. Refer to the appropriate SOP windows - **Sample** and **Flow settings** - for more details.

The **Edit settings** button will display the standard **Edit result** window where the measurement details and parameters can be inspected and edited.

③ Traces tab

The **Traces** tab displays the traces available for viewing on the plot. The traces available will be determined by the flow-mode system arrangement and the measurement signals analysed.

The available traces could be the **Z-average mean, Intensity, RI** (Refractive Index) and **UV**. When **RI** is selected using the drop-down menu, the **Delay Volume, Offset** and **Gain** values will be displayed; These are the values entered in the **External Inputs** SOP window when the SOP was originally performed. Refer to the **Flow-mode SOPs** section.

Note that the default **External Inputs** values are defined in the **Options** window, explained later in this chapter.

④ Peaks tab

The **Peaks** tab displays all the details about the calculated peaks shown in the flow-mode plot. To display a peak's details either select a blue peak bar in the plot, or select one of the peaks identified at the top of the peak tab; the details in the tab will reflect the data appropriate to that peak. The selected peak will be highlighted on the plot. (The details displayed are the same as shown on the **Chromatogram summary** report).

The features of the **Peaks** tab are:

Peaks identified box	This shows all the peaks identified in the measurement and plot. Selecting a peak will display its details in the accompanying boxes. To show or hide a peak in the plot, select or de-select the check box alongside the peak.
Hide all	Press this button to hide all identified peaks. Peaks can be displayed again by selecting the peak check box.
Remove	This will remove the highlighted peak in the Peaks identified box. When pressed, a message is displayed to warn that the peak will be removed permanently.
Start vol. / End vol.	This shows, in mL, the start and end volume points for the peak.
Peak analysis details	Displays the results from the measurement: Molecular weight, Z-average, Width and Intensity area .

Controlling the flow-mode measurement via an SOP

Once the Flow-mode option has been connected to the required size exclusion chromatography (**SEC**) system, a measurement can be made. This is done using an SOP, or a manual measurement in the usual manner.

The **Flow**-mode SOP enables measurements to be performed on a flowing sample stream. The scattered light intensity and hydrodynamic diameter can be plotted as a trend, and, with the addition of optional external input windows, information from external measuring sources, such as a refractive index detector, can be monitored and displayed.

An SOP can be configured to control all settings for the accessories automatically.

A **Flow-mode** measurement follows the same SOP format as performed when doing a normal **size** measurement, with a few exceptions. When a size measurement is chosen some extra dialogues - **External Input 1/2** and **Flow settings** - will be included in the SOP selections.

The **SOP Editor** and setup is described in full in the **Zetasizer Nano user manual**. Most of the SOP sections are common to Measurement types, and these are described in the above manual. The other SOP sections are specific to the flow measurement being performed; these are described below. Also note that some of the other dialogue pages will alter slightly to accommodate extra parameters necessary to perform the flow measurement.

Creating or editing an SOP - Measurement type selection

- To **create** a new SOP, select **File-New SOP**.
To **edit** an existing SOP, choose **Open-SOP**.
This will open up the SOP editor. The SOP editor consists of several dialogues that can be stepped through by using the **Next** arrow button.
- Complete the SOP editor as described in the Zetasizer Nano user manual.
- Once the SOP has been created, press **Finish** and save the new SOP.




Note

To improve the stability of temperature control when using the system as a Gel Permeation Chromatography (GPC) detector, the measurement temperature should be set to be 5°C or greater, above or below the ambient temperature.

Measurement type options

Select a **Flow** measurement type then complete the SOP creation as required.



External input 1 and 2



Measurement parameters can be set for each **external input** used. The same parameter options are available for each input.

Settings	Description
	With the Enable input check box selected, the input conditions required for the external signal can be setup.
Input name	Use Input name to name the type of input. Use a name that represents the signal or reading being inputted. A name is mandatory once the check box has been selected, otherwise another SOP window cannot be selected.
Input type / Path length	Input type (the kind of detector the instrument is connected to) and Path length (mm) are options that only require selection when using the Research features .

Settings	Description
Offset / Gain	The Offset and Gain are mathematical parameters (values) that are needed to convert the incoming signal (in volts from the external measurement device) to the measurement parameter required.




ill 7696

Consult the external detector documentation to ascertain the **Offset** and **Gain** values required and input these into their respective entry boxes.

Delay volume	<p>The Delay volume is the volume of liquid contained in the cells and connecting tubing between the external detector and the Zetasizer Nano. This parameter ensures that data recorded by the Zetasizer can be over-plotted from the same elution point.</p> <p>This value will be obtained from the external detecting device (i.e. chromatographic column) that is used in the flow-mode measurement.</p> <p>The default settings, for name, delay volume, offset and gain, can be adjusted using the Tools-Options-External Inputs window.</p>
--------------	--

Measurement



Settings

Measurement duration


Description

The measurement duration setting may affect the accuracy and repeatability of the results.

In **Measurement duration**, input the total measurement time or volume amount required, and adjust the units to suit: **Time** or **Volume**.

The **Run duration** value determines the length of each individual measurement within the experiment.

Flow settings



Settings	Description
Flow rate	Input the flow rate of the sample through the instrument and connecting tubing. This value is taken from the external detecting device (i.e. chromatographic column) that is used in the flow-mode measurement.

For all other SOP windows, refer to the size SOP section.


The Flow-mode measurement display

The Flow-mode measurement displays are virtually identical to those shown when performing a standard size measurement. The only difference being the inclusion of different view options on the Result tab, i.e. **Flow trace vs Volume** as shown below.

The standard tabs are explained in the Zetasizer Nano manual.

Result tab (1st tab)

The result tab will show the result obtained as the measurement progresses. The result view will be updated after every run of the measurement.



The Result tab is named after the result view chosen, the result view shown above is **Flow trace vs Time**. Different views can be selected by right-clicking on the graph and selecting from the list displayed.

The views available are: **Count Rate**, **Correlation Function**, **Flow trace vs Volume** and **Flow trace vs Time** and **Monitor**.

■ **Monitor** enables the count rate signals to viewed before a measurement is run.


Note that the monitor only displays information before the measurement has actually been started. When the measurement is started, the flow trace views should be used.

Displaying the flow-mode report

To display a titration report, select a 'Flow' type measurement record and then select the appropriate **report tab**. The report will show all appropriate measurement information for that record.

Standard report - flow-mode measurements

The standard Flow-mode report, **Chromatogram Summary (M)**, gives the same information as seen in a standard size report, plus additional information relating to the flow duration and rate used.



The three main peaks in the measurement will also be shown, displaying the sample intensity, width, molecular weight and start and end flow volumes.

Additionally a the flow trace result graph will be displayed. This can be viewed either in **Time** or **Volume**.

Microrheology

Introduction

This chapter gives an overview of the **microrheology** option in the Zetasizer software. It can be used with either a Zetasizer Nano ZS or ZSP. The **microrheology** option allows the measurement of the viscoelastic modulus of samples within the linear viscoelastic region.

ZEN5600 Zetasizer Nano ZSP instrument

ZEN3600 Zetasizer Nano ZS instrument

This chapter provides the following information:

- About DLS microrheology
- Microrheology measurement process
- Microrheology SOPs
- Sample preparation advice
- Running microrheology measurements
- Analyzing microrheology data

About DLS microrheology

The term for the measurement protocol that a Microrheology optioned Zetasizer Nano instrument performs is DLS Microrheology (μ Rh). This technique uses dynamic light scattering to determine the rheological properties of low viscosity and weakly structured samples by measuring the motion of an embedded tracer particle.

DLS Microrheology uses tracer probe particles to measure the relationship between stress and deformation. Both DLS and ELS (Electrophoretic Light Scattering) measurements are utilized in the microrheology application as part of measurement and optimization protocol.

Advantages of DLS Microrheology (μ Rh) are:

- Rheological characterization of low viscosity, weakly-structured samples. Via access to very high frequency (short time) dynamics which is highly relevant for dilute samples or weak structures.
- Very small sample volumes can be used (approx. $12\mu\text{L}$). This is particularly suitable for precious protein samples.
- Provides rheological parameters - G' , G'' , η^* . Data can be verified using same sample measured on a rotational rheometer where measurement ranges overlap. Data can be exported and used in the Malvern rSpace software.
- Extends viscoelastic measurements into ranges inaccessible by mechanical rheometry techniques.
- Fast measurements are possible, with all frequencies effectively sampled simultaneously.
- Applications for viscosity of protein solutions and onset of protein aggregation.

Microrheology measurement process

With the cell loaded into the instrument, and the SOP configured, a measurement can be performed.


To make a microrheology measurement, only one measurement is technically needed - a size measurement used by the software to calculate the relaxation times, and therefore the viscosity and moduli values.

However, in order to make reliable measurements, the full process consists of making three distinct measurements, in this sequence:

1. **Tracer compatibility measurement (zeta potential)**
2. **Tracer concentration measurement (size)**
3. **Microrheology measurement**

Steps 1 and 2 are optimization phases. These are zeta potential and size measurements on the tracer particle and subsequently the sample to ensure all requirements are met before performing the microrheology measurement itself. These stages are not always necessary and depend upon the measurement setup and data available. It is possible to proceed directly to the microrheology measurement (3) if required.

The optimization stages are shown in the flow charts below:



Microrheology SOPs

Creating a microrheology SOP or manual measurement is very similar to that of a standard SOP, except that one of the **Microrheology** measurement types is selected:

Measurement type	Description
Microrheology optimization – Tracer compatibility	<p>Zeta Potential measurement to check there is no significant interaction between the tracer particles and sample matrix.</p> <p>The tracer is initially measured in the sample buffer/dispersant (or solvent) alone with no sample added, sample is then added and a second measurement performed. The two zeta potential results are then compared - a small relative change in zeta potential in the presence of the sample indicates minimal interaction between the sample and tracer particles.</p>
Microrheology optimization – Tracer concentration	<p>Size measurement to evaluate the optimum tracer concentration.</p> <p>Used to assess the relative intensity of tracer particle scattering against sample matrix scattering. This ratio is then used to predict the tracer concentration required to give dominant scattering (95% of total scattering intensity).</p>
Microrheology measurement	<p>Measure Correlation function of tracer particles and extract Microrheology data from Mean Square Displacement (MSD) plot.</p>



Note:

For a description of the other measurement parameters, refer to the help system.

Tracer settings

Tracer SOP settings are configurable within the optimization measurement types:

Option	Description
Tracer details	The parameters of the tracer that will be used in the measurement.
Select	Highlight the tracer required and click OK . The selected measurement will use the properties of that tracer in the test. Only a microrheology measurement requires the nominal tracer size, but each optimization test should ideally be linked to the appropriate tracer type and chemistry. A new tracer entry can be created by clicking on an empty row and completing the required details.
Edit	Enables the nominal tracer size to be edited. For Malvern supplied tracer particles this is usually the nominal diameter given on the bottle but ideally the hydrodynamic size of the tracer in the sample solvent should be used. This can be determined from a standard size measurement.
Notes	Add any additional information as required.
Tracer concentration (µl/ml)	Set the concentration value of the tracer (available in tracer concentration measurements only).

Advanced settings

Several Advanced SOP settings are available for microrheology - the option shown depends on the measurement type selected:

Measurement type	Option	Description
Tracer Compatibility	Acceptable Zeta Difference (mV)	This optimization test compares the two measured zeta potential results and tests to see if sample measurement is within the set tolerance (mV) of the tracer measurement. The default value is 5mV.
Tracer Concentration	Extend duration for high viscosities	Enable this option to automatically extend the measurement time when particles with high viscosities are present. Useful when measuring the forward angle.
Microrheology	MSD smoothing window size	The (odd) number of points in the smoothing window used to smooth the MSD (Mean Squared Displacement) data. The available settings are: No smoothing, 3, 5, 7, 9.

Sample preparation advice

Correct preparation of the tracer and sample elements is important in ensuring the reliability and repeatability of the results. The concentration of tracer samples depends to a great extent on the concentration of your sample in the dispersant.

The optimization steps for **Tracer Compatibility** and **Tracer Concentration** are useful tools for assessing whether the tracer chemistry is likely to be suitable for a particular sample, and if so the minimum tracer concentration that can be used.

**Note**

If this information is already known then the user can proceed straight to the **microrheology measurement**.

Following are some important considerations for the preparation.

Tracer compatibility

For the initial tracer compatibility tests (zeta potential), the following starting concentrations are suggested.

- **Tracer particles in dispersant** - 5 μ l neat tracer particle suspension to 10ml dispersant (solvent) (0.5 μ l/ml).
- **Tracer particles in presence of sample** - add one drop of the dispersed (diluted) sample (approx. 0.05ml) to approximately 5ml of the previously prepared tracer dispersion allowing sufficient mixing and equilibration time prior to measurement.

Acceptable Zeta Difference

Testing the difference in the zeta potential between the tracer and measurement sample is important because large changes would indicate that the tracer particle and the sample are interacting. With proteins, that could mean that the proteins are adhering to the tracer's surface and therefore not only changing the tracer's zeta potential but also the size of the tracer. The tracer then becomes part of the system and directly influences the rheology.

Therefore to avoid this interaction the amount of sample added with respect to the tracer must be controlled. This is done by measuring the **sample** zeta potential measurement against the **tracer** zeta potential. The resultant measurement and comparison must show that the **sample** measurement is within a set tolerance of **tracer** measurement.

This comparison is termed the **Acceptable Zeta Difference**.

For the microrheology measurement to measure reliably the **Acceptable Zeta Difference** is **5mV or less**.

**Note**

The complete dispersion of the tracer in the sample can take a significant amount of time. For certain systems, gentle mixing (for example using a sample roller) will help to disperse the tracer particles. However, surface interactions can take several hours to manifest themselves.

During the measurement this difference is checked and the software advises the user.

There are two possible causes for the zeta potential showing significant differences between the tracer and the sample / tracer systems.

- **The tracer concentration is not high enough in the sample**

If the tracer concentration is too low, the zeta potential result will be dominated by the sample, not the tracer, adding more tracer particles and re-testing can show if the difference is due to concentration rather than interaction.

- **The tracer and sample are interacting**

If the tracer and sample are interacting, adding more tracer will not improve the difference in zeta potential (unless so much is added that it completely dominates the zeta potential measurement). Another tracer should be tested to see if it is more suitable for the application.

Tracer concentration

Tracer concentration Size measurements are performed to check that the tracer particles' scattering completely hides the sample scattering; therefore the size of the tracer particle should be larger than the sample.

- **Suggested concentration of tracer particles**

Add approximately 5 μ l of neat tracer particle suspension to 1ml of sample (or a relative amount depending on the type of cell being used) and ensure the particles are well dispersed. The exact concentration should then be entered in the **Tracer Concentration** field in the SOP's Tracer properties section.

- **Filtering**

It is suggested that to get the most reproducible result, the tracers are filtered with appropriately sized (i.e. filter size larger than the tracers) filter.

Running microrheology measurements

Use either a manual measurement or SOP to run one of the microrheology measurement types.

- **Tracer compatibility measurement (zeta potential)**
- **Tracer concentration measurement (size)**
- **Microrheology measurement**

During the measurement's progress the display is similar to that of a standard zeta potential or size measurement, with the addition of extra information pop-up windows. This section provides a run through of each measurement type.

Tracer compatibility measurement

The **tracer compatibility measurement** type is used to assess the suitability of tracer particle chemistry in order to minimize particle-matrix interactions by measuring:

1. The zeta potential of the tracer particles in the dispersant (solvent) being used for the Microrheology test.
2. The zeta potential of the tracer particles in the dispersant AND sample being used for the Microrheology test.

The test compares the two measured zeta potential results and tests to see if sample measurement is within a set tolerance of the tracer measurement (the 'Acceptable Zeta Difference (mV)'). The software default value for this difference is 5mV, which can be changed in the advanced settings.

The zeta potential of a colloidal particle will change if the surface properties of the particle change e.g. if components of the dispersant matrix, such as polymer or protein molecules, adsorbed onto the particle surface. If the zeta potential of the tracer in the presence of the sample is significantly different from the value obtained in the dispersant (solvent), this may indicate strong particle-matrix interactions, that will likely affect the rheological data obtained.


Running the measurement

This section assumes you have already set up the SOP or manual measurement parameters appropriately.


**Note**

Only use a zeta potential cell for this procedure: **ZEN1002 Universal dip cell** or **DTS1070 Disposable folded capillary cell**.


1. Fill a zeta potential cell with the dispersant/buffer and add a diluted suspension of the tracer particles.
2. Start the measurement - an initial measurement of the zeta potential of the tracer is performed and the following message is displayed:



3. Click **Next**. The **tracer particle zeta potential** measurement is made and the result shown:



4. Follow the instructions given and click **Next**.
5. The variance from the **Acceptable Zeta Difference** is displayed:



- If the difference is **less than 5mV** then the measurement has been successful - proceed to make a **Tracer concentration measurement**.
- If the difference is **more than 5mV** the software will provide sample preparation advice and advise on repeating the measurement.

Tracer concentration measurement

The **tracer concentration measurement** assesses the relative intensity of tracer particle scattering against sample matrix scattering. This ratio is then used to predict the tracer concentration required to give dominant scattering (95% of total scattering intensity).

The measurement utilizes a standard size measurement and records the scattering intensity from the tracer particles relative to that of the sample. Using the relative scattering ratio and concentration entered, the software will then calculate the concentration of tracer particles required to give a relative scattering intensity of 95% (19:1 intensity ratio), which is the minimum recommended to ensure that tracer scattering will dominate the measurement.


This calculation is made on the basis that for single scattering there should be a linear correlation between concentration and intensity.

Running the measurement


This section assumes you have already set up the SOP or manual measurement parameters appropriately.

1. Add approximately 5 μ l of neat tracer particle suspension to 1ml of sample (or a relative amount depending on the type of cell being used) and ensure the particles are well dispersed. Enter the exact concentration in the **Sample-Tracer** window under **Measurement Parameters**.

2. Once prepared transfer the sample to a size capable measurement cell and start the measurement. The following information is shown:



3. Click **Next** - the system performs a size measurement and then provides a summary of the result, showing the recommended concentration value:



Microrheology measurement

The Microrheology measurement is a single size measurement from which the MSD and rheological parameters are calculated.

The process measures the correlation function of the probe (tracer) particles in the sample matrix under test. The Mean Square Displacement (MSD) of probe particles is then calculated, from which the rheological properties of the sample are subsequently evaluated (using the Generalized Stokes-Einstein Relation).

The only extra information the user requires for the microrheology test is the **hydrodynamic size** of the tracer particle in their sample solvent or buffer. This step is not included as an optimization step in the microrheology application as it can be determined from a standard size measurement of the tracer in the solvent.

Some tracer particles are already set-up in the Sample-Tracer window under Tracer Details. These include tracer particles supplied by Malvern and can be 'edited' or 'added to' by the user as required. For Malvern supplied particles the Nominal Tracer Size is usually the size on the bottle which may differ from the hydrodynamic size in the sample solvent. This value can be modified by ticking the Edit box. It is this value which is used in the microrheology calculation.

An optional smoothing window can be applied to the Mean Square Displacement prior to conversion in to rheological parameters in order to improve the quality of the final data. This can be found in the **Measurement-Advanced** window and is set-up to use a 5 point smoothing window as default.

**Note:**

Both the **Nominal Tracer Size** and **Smoothing Window** can be edited after the test if required with a new result entry created with the modified values.

Analyzing microrheology data

After running a microrheology measurement, analyse the data in any of the ways discussed in this section:


- The microrheology reports
- Microrheology utilities
- Exporting the data for analysis in another software tool, such as rSpace

Microrheology reports

There are a number of standard reports available for reviewing the measurement results:

- **Microrheology:** Correlogram and Mean Squared Displacement (MSD).
- **Rheological properties:** Complex viscosity (η^*), Viscoelastic Modulus (G), Creep compliance (J).

To display a microrheology report, select a **Microrheology measurement** record and then select the appropriate **report tab**. The report will show all appropriate measurement information for that record. This is an example of the Mean Squared Displacement Report (MSD):



Mean Squared Displacement (MSD): a representation of the movement of the tracer particles within the sample. A purely viscous (Newtonian) sample will show a straight line.

Correlogram: shows the base data from which the microrheology result is calculated. The correlation function displayed is a measurement of the amount of movement the probe particles make within the test sample over a range of time scales.

Complex viscosity (η^*): shows the relationship of the complex viscosity to the shear viscosity using the Cox-Merz rule.

This rule is generally only fully applicable for simple systems, and that the differences between complex viscosity and shear viscosity increases as the sample structure becomes more complex.

Viscoelastic Modulus (G): The Moduli tab shows the viscoelasticity of the sample, it shows two sets of data on the chart - G' , the elastic (storage) component and G'' , the viscous (loss) component.


In a microrheology measurement it is expected that the viscous component will be dominant for at least most of the measured frequency range. This is because the technique relies on the tracer particles being able to move and therefore produce a correlogram.

For example; a gel system where the elastic component (G') is dominant for all of the frequency range, will hold the tracer particles still within the gel matrix, which means that the correlogram will be close to a flat line, and the MSD will be effectively 0.

Creep compliance (J): Creep compliance is another way of viewing the viscoelasticity of a sample.

Microrheology Utilities

The Microrheology Utilities allow you to analyze data from a microrheology measurement. Right click a Microrheology record and click **Microrheology Utilities** or choose **Tools-Utilities-Microrheology Utilities**.



The following tab selections are available:

- Mean square displacement
- Viscoelastic Moduli
- Complex viscosity

Each tab page contains controls that allow different models to be fitted to the microrheology data. It is important to note however that not all models will be applicable to all the datasets. Due to the measurement inconsistencies at very short timescales, there will be areas of the data that will not be appropriate to use with the models. In order to change the area used, and displayed in the Fit results area, left-click and drag the range pointers (the red triangles on the x axis) to an appropriate point on the chart. To export the microrheology data, refer to the Microrheology data topic.

Exporting microrheology data

Following the completion of the measurement, the rheology data can be exported from the Zetasizer Nano software as an .xml or a .csv file.

**Note**

The .csv format should be used if the results are to be imported into the Malvern rSpace software.

Exporting the microrheology data (.xml/.csv) - using the Microrheology utilities:

1. Select a microrheology measurement record.
2. Select **Tools-Utilities- Microrheology utilities**.
3. The Microrheology data processing window is displayed. Choose **Save** and select the export destination. Alter the file extension to .xml file or .csv as required.
4. Click **OK** to export the data. The parameters exported are:

Sample Name	Date
File Name	Lag times (μ s)
Times (μ s)	Angular Frequency (rad/s)
Creep Compliance: Plot against Times	Mean Square Displacement: Plot against Times
Channel Values: The Correlogram, plot against Lag times	Complex Viscosity Plot against angular frequency
Shear modulus (Viscous component): Plot against angular frequency	Shear modulus (Elastic component): Elastic component, plot against angular frequency

**Note**

This data can be imported into rSpace software using an appropriate sequence.

Advanced protein features

Introduction

This chapter gives an overview of the **Advanced protein features** option in the Zetasizer software.

- With the Zetasizer Nano ZSP it is now possible to achieve the best possible measurement of **protein mobilities**. This is achieved with the combination of the following features:
- A system with sufficient sensitivity to measure the low count rates and low electrophoretic mobilities associated with protein samples.
- A measurement technique that minimises the risk of protein aggregation and minimise the amount of material required. This is the **Diffusion barrier** measurement technique; a description of this follows.
- A measurement process that reduces the risk of aggregation but also capable of identifying any aggregates that are formed and assesses the quality of the measurement.

Measurement process


The protein mobility measurement combines **size** and **zeta** potential measurements to check that no protein aggregates are forming during the measurement.

The mobility measurements will be performed in **groups** of sub-runs. This is to allow for periods of cooling to be applied at stages within the measurement and so reduce any chances of aggregation happening through over-heating of the sample.

To further reduce any risk of aggregation, an automated measurement optimization will check and reduce the voltage required. This is done by monitoring the sample conductivity.

The basic steps involved in a measurement are:

- The first step is the **thermal equilibration delay**, performed in order for the sample and cell to properly equilibrate with the Zetasizer cell holder.
- An optional **pre-mobility size** measurement is completed so that the user can check that the sample is not aggregated.
This result is recorded into the protein mobility parent record as a child record of type **size**.
- A **zeta potential** measurement follows for measuring the mobility.
A normal zeta potential measurement consists of a number of sub-runs, at the end of which the final result (the average over all recorded sub runs) is reported into the record view.
- During the measurement proper, running a large number of sub-runs sequentially will significantly increase the risk of Joule heating of the sample. to prevent this the protein mobility measurement is split into the smaller **groups** of sub-runs with a delay added between each group to allow the sample to relax.



- Once the measurement starts the **Expert** system will check and warn if the sample has will be subject to levels of Joule heating or Polarisation/field estimation issues that might affect the measurement accuracy. Appropriate advice is given, for such cases, on the **Expert advice tab** on the **live display** and after the measurement on the record view. iii iii8801


- Once all of the groups have been recorded the **post-mobility size** measurement is performed in order to determine the state of the sample after the electrophoresis measurement.
- This completes the measurement process.

Diffusion barrier measurement technique

Protein mobility measurements should be used in combination with the **diffusion barrier** measurement technique to further protect the sample.

It has been found that most of the damage that happens to protein samples, during mobility measurements, happens at the electrodes. The diffusion barrier technique protects the protein from damage by introducing a physical distance between the sample and the electrodes. This is done by holding a plug of measurement sample within a larger volume of the same buffer that contains no protein. This separates the protein from the electrodes and prevents any protein-electrode interaction and sample damage.

- Ⓐ. Standard measurement volume
Measurement sample covers electrodes
- Ⓑ. Diffusion barrier measurement volume
Sample plug shown with buffer between plug and electrodes



iii iii8802

The sample sits at the bottom of the cell if the cell is handled carefully.

Advantages of this technique are that:

- Many more measurements can therefore be made before aggregates start to appear.
- Only small volumes (20-50 μ l) are required.

Controlling a protein mobility measurement via an SOP

A **protein mobility** measurement is divided into two separate measurement parts - **zeta potential** and **size**. The **protein mobility** measurement follows the same SOP format as performed when doing a normal zeta potential or size measurement, with a few exceptions.

The **SOP editor** and setup is described in full in the **Zetasizer Nano user manual**. Most of the SOP sections are common to Measurement types, and these are described in the above manual. The other SOP sections are specific to the protein mobility measurement being performed; these are described below. Also note that some of the other dialogue pages will alter slightly to accommodate extra parameters necessary to perform the microrheology measurement.

Creating or editing an SOP - Measurement type selection


- To create a new SOP, select **File-New SOP**. This will open up the SOP editor. The SOP editor consists of several dialogues that can be stepped through by using the **Next** arrow button.
(To edit an existing SOP, choose **Open-SOP** instead.)
- Complete the SOP editor as described in the Zetasizer Nano user manual.
- Once the SOP has been created, press **Finish** and save the new SOP.

The various SOP dialogues are described below.

Measurement type options

Select a **protein mobility** measurement type then complete the SOP creation as required.


Measurement



Settings	Description
Sub-runs per group	<p>With Automatic selected, the software will automatically determine the number of sub-runs required per measurement. This will be suitable for the majority of samples and can simply be left as the default. With manual selected the measurement will use the user defined Number of runs setting.</p> <p>The time may be reduced for the measurement of a latex standard, or increased to improve the repeatability of the measurement of particularly polydisperse samples. All the individual runs are accumulated and then averaged to give a final zeta potential result. Therefore the more runs performed the better the repeatability. Naturally the more runs selected the longer the duration of the complete measurement.</p>
Groups per measurement	<p>Measurements of protein mobility are run as groups of a limited number of sub-runs. The results from the groups are combined in the final result. This option defines the total number of groups that are combined to make up the result. Add a delay between the measurements of each group in the Delay between groups entry box if required. A typical delay between groups would be 180s.</p>

Running a large number of sub-runs sequentially significantly increases the risk of Joule heating of the sample, so the protein mobility measurement is split into the smaller **groups** of sub-runs with a delay between groups to allow the sample to relax.

Size measurement



Full details of the size measurement window will be found in the main Zetasizer Nano user manual. A brief overview follows.

Settings	Description
Perform a size measurement...	With the Perform a size measurement before and after the Protein mobility measurement check box selected a size measurement will also be performed as well as the standard zeta potential measurement.
Angle of detection	The measurement angle is fixed to 13 degree forward scatter. This is to ensure that the size is measured at the same angle as the zeta potential.
Measurement duration	The Measurement duration setting can affect the accuracy and repeatability of the size results.
Measurements	This option allows a sample to be measured more than once; to investigate the effect on particle size over time, or to prove repeatability.

Settings	Description
Partial results	If it is likely that a measurement will not produce a correlation function that can be analysed, then the data collected can still be saved by selecting the Allow results to be saved containing correlation data only check box.

For all other protein mobility SOP windows, refer to the size SOP section in the main user manual.

Performing the measurement

The protein mobility measurement follows the same basic measurement sequence as a standard zeta potential measurement.

The measurement will be performed in **groups**, as defined in the SOP. These are groups of a few sub-runs with a delay between each group to allow the sample to relax.

The protein mobility measurement displays are identical to those shown when performing a standard zeta potential measurement.

The measurement and a description of the measurement tabs is explained in the Zetasizer Nano user manual.

During the measurement the data is stored as a **'parent'** protein mobility zeta potential record, with **'child'** records relating to the individual zeta potential and size measurements made during each measurement group.

Record	Type	Measurement Date and Time	Sample Name	T °C	Z-Ave r.nm	Pdl	Mean Count F kcps
50	Protein Mobility	01 August 2012 15:19:31		1			
54	Protein Mobility	01 August 2012 15:27:18		1			
55	Protein Mobility	01 August 2012 15:27:37		1			
56	Size	01 August 2012 15:27:59	1 1	25.0	50.00	0.002	
57	Size	01 August 2012 15:28:02	1 2	25.0	50.00	0.002	
58	Size	01 August 2012 15:28:05	1 3	25.0	50.00	0.002	
59	Zeta	01 August 2012 15:28:06		1 25.0			
60	Size	01 August 2012 15:32:34	1 1	25.0	50.00	0.002	
61	Size	01 August 2012 15:32:37	1 2	25.0	50.00	0.002	
62	Size	01 August 2012 15:32:40	1 3	25.0	50.00	0.002	
Mean 55							
Std Dev							

Interpreting the results

The record view

As part of each group measurement, the count rate is measured and a rolling average taken as the measurement proceeds.

Aggregates are, generally, characterised by a much larger particle size than the native protein and if any aggregates are present then large changes in the measured count rate in each group will be observed. This will happen even if there are only tiny fractions of aggregates present with respect to the fraction of protein.

Therefore once all of the measurement groups have been recorded, any groups that have a significantly higher count rate than the average will be removed from the averaging procedure.

The mobility over the remaining groups is then calculated and reported into the record view as a **zeta** potential child record of the **protein mobility** parent record.



Note

The **Expert** system will warn the user if **less than ten** groups are left once aggregated groups have been removed.

The average over un-aggregated groups will yield the mobility of only the protein itself and not its aggregates, but the number of groups included in the measurement will aid the user in assessing whether the measurement of further aliquots is required.

- The **pre-mobility size** records tell the user whether the sample is aggregated prior to the mobility measurement.
If it is then the aggregates will, typically, scatter far more light than the native protein and the aggregate mobility will dominate the result.
- The **zeta** potential record is the measurement of the mobility of the protein with the aggregated groups removed prior to calculating the final result and the **Expert** system will give advice on the quality of this measurement with regards to field effects such as Joule heating.
- Finally the **post-mobility size** measurement characterises the sample after the mobility measurement.

This is especially important when the Diffusion Barrier is used, where the sample may not have reached the electrodes; subsequently no electrode aggregation should have occurred and the sample can be retrieved for further analysis.

Displaying the protein mobility measurement report

To display a protein mobility report, select a **Protein mobility group** type measurement record and then select the appropriate **report tab**. The report will show all appropriate measurement information for that record.

Standard report - protein mobility measurements

Once a protein mobility measurement is completed there are a number of standard reports available for reviewing the measurement results.

Protein Mobility (M)	
Sample Name: 13	
SOP Name: mansettings.nano	
File Name: MicroRheoData.dts	Dispersant Name: Water
Record Number: 62	Dispersant RI: 1.330
Date and Time: 01 August 2012 15:32:40	Viscosity (cP): 0.8872
	Dispersant Dielectric Constant: 78.5
Temperature (°C): 25.0	Total Sub-Runs:
Count Rate (kcps):	Sub-run groups: 0
Cell Description: Disposable sizing cuvette	Measurement Position (mm): 4.65
Attenuator: 11	
Result quality : Not a zeta measurement	
Mobility ($\mu\text{mcm/Vs}$):	
Mobility Deviation ($\mu\text{mcm/Vs}$):	
Conductivity (mS/cm):	

The main reports are **Protein mobility** and **Protein analysis**. An example of the Protein mobility report is shown above. The report gives similar information as seen in a standard size and zeta potential report, plus additional information relating to the protein mobility measurement itself.


Calculators tool

One of the three basic functions of the Zetasizer Nano series of instruments is its ability to perform accurate measurement of a sample's molecular weight. By measuring the sample scattering intensity over a range of concentrations and entering the necessary sample parameters, the molecular weight can be determined.

If the hydrodynamic diameter is also measured from one of these concentrations, the molecule shape or conformation can also be estimated.

The **Calculators** tool enables the calculation of the molecular weight and also provides other calculation tools.

To access the protein utilities select **Tools-Calculators**. The following tab selections are available:



- **SLS Debye plot**

The protein utilities function enables a static light scattering **Debye plot** to be constructed from freely available information, or from a record generated from actual data.

- **Molecular weight (MW) and shape estimates**

A 'what if' calculation can be performed. If any two parameters from the **molecular weight**, hydrodynamic size and conformation (**shape**) are known, then the third - either a Shape, Hydrodynamic diameter or molecular weight parameter - can be estimated.

- **Concentration utilities**

This window contains features to establish the concentration and scattering levels that are expected to give the input parameters.

- **Scattering functions**

A plot can be generated by entering the measurement data.

- **Protein charge & $f(Ka)$**

A tool to calculate $f(Ka)$ from the Henry equation using known size and ionic strength. This tool uses the Ohshima equation for monovalent salts. The

second calculation calculates protein charge from the electrophoretic mobility and Stokes radius.

■ **Virial Diameter**

A tool to calculate the virial diameter, also called the thermodynamic diameter, from the measured molecular weight and 2nd virial coefficient (A₂).

■ **DLS Debye plot**

A 'dynamic Debye plot' of measured hydrodynamic size as a function of sample concentration to calculate the true hydrodynamic radius and the dynamic virial coefficient.

■ **Interparticle distance**

A tool to calculate the distance between the particles based on their concentration and molecular weight. Also estimates the thickness of the electrostatic layer based on protein charge and ionic strength.

■ **Mixture Viscosity**

A simple tool that calculates the overall viscosity of a mixture of solutions based on a volume weighted mean viscosity.


■ **Oligomer ratios**

A tool that estimates the ratio between monomer and dimer in a peak that contains both where the size of each is known. This is based on work published by Malvern Instruments entitled "Dynamic light scattering as a relative tool for assessing the molecular integrity and stability of monoclonal antibodies" by Nobbmann U et al. The reference is published in Biotechnology and Genetic engineering Reviews, 2007, Vol 24 pp117-128.

In each case the effect of changing any of the input parameters can be seen instantly in all of the derived parameters.

SLS Debye plot

As mentioned in the introduction, the **Calculators** tool includes the ability to generate a Debye plot, using inputted rather than measured data.



This feature can be useful for various reasons, for example:

- By combining individual measurements, one single Debye plot can be generated.
- A Debye plot can be created by entering the concentration points from an existing measurement, then adding additional concentration points.
- Any of the parameters in the window box can be changed; the other parameters will be instantly recalculated. This can be used to investigate the sensitivity of the result to changes in any parameter.

For example, by first entering the concentration points from an existing measurement, a sample parameter - e.g. sample temperature - can be altered and the effect immediately observed on the Debye plot. This saves time in performing the original measurement again at the different temperature.

The format of the plot can be altered by changing the drop-down menu in the top right hand corner of the window to **Chart properties**. The individual properties can then be altered in the table below.

**Note**

The Debye plot uses the reduced Rayleigh scattering equation.

Adding & editing sample parameters & table data

To generate an SLS Debye plot, the sample parameters and table concentrations have to be entered.

To access the **Debye plot** select **Tools-Calculators** and then the **SLS Debye plot** tab. Select the **Table** tab to begin entering data into the table.

Concentration		Intensity	Residual Intensity	Ro
(mg/mL)	(g/mL)	(kcps)	(kcps)	(1/cm)
0.000	0.000000		0.0	-1.36685e-
10.000	0.010000		10.0	-1.28747e-
20.000	0.020000		20.0	-1.20810e-
30.000	0.030000		30.0	-1.12872e-
40.000	0.040000		40.0	-1.04935e-

- To **define a new** concentration, press the **Add...** button. The table **Input values** window will be displayed.


Specify the **Concentration** and **Sample intensity** values - either new values or ones taken from an existing measurement.

Input Values

Concentration (mg/ml):

Sample Intensity (kcps):

- To **modify** a concentration, select it from the list and press the **Modify...** button. The **Input values** window is displayed, allowing the parameters to be changed.
- A concentration can be **deleted** by selecting the concentration from the list and pressing the **Delete...** button.
- Select the **Graph** tab to see the resultant **Debye plot**.



- The table values and graph plot can be subsequently altered by changing the **Sample, Data** and **System** parameters in the measurement parameters table on the right of the window. These parameters are described in the following section.

Measurement parameters table

Once all the concentration values have been added into the table, the measurement parameters table can be used to alter the result and Debye plot.

To view and alter each parameter setting, click on the plus \blacktriangle sign next to each parameter group to open the list (Click on the minus \blacktriangleright symbol to close the list).

Default parameters are in “normal” type, altered parameters will be made **bold**.

Debye Calculator	
\blacktriangle Data	
\blacktriangleright Data points	5
\blacktriangle Sample	
Dark Count	0
dn/dc (ml/g)	0.105
Notes	
Solvent Ref. index	1.499
Solvent count rate	172.2
ShapeModel	Small Molecule (no shape correc
\blacktriangle System	
Wavelength (nm)	633
Scattering angle	175
Toluene count rate	172.2
Dark Count	
The value to use for correcting the scattering using Dark Count	

The parameters are described below.

Data

This indicates the data entries that have been entered into the Debye plot table, see below.

Sample

- **Sample dn/dc (ml/g)**

This is the specific refractive index increment; the change in refractive index as a function of the change in concentration.

- **Experiment notes**

Used to record specific details about the experiment or calculation performed.

- **Hydrodynamic radius (nm)**

The radius as measured using dynamic light scattering.

- **Solvent Ref. Index (Refractive index)**

The refractive index of the solvent used.

- **Solvent count rate (kcounts)**

The count rate used to calculate the Rayleigh ratio (R_0) is the 'residual' count rate, which is derived by subtracting the solvent count rate from the sample count rate.

- **Shape Model**

The shape model that is used to estimate the radius of gyration from the hydrodynamic radius, and therefore calculate the angle dependent effects on KC/RoP for particles of sizes outside the Rayleigh region (Diameter > ~50nm).

System

- **Wavelength (nm)**

The wavelength of the laser used in the Zetasizer Nano instruments, or for the measurement. Either 632.8nm 'red' or 532nm 'green' laser wavelengths are available.

- **Scattering angle (degrees)**

The scattering angle of the optics unit - either 173°, or 90° for the 'classical' optics arrangement.

- **Toluene count rate (kcounts)**

The scattering count rate of the toluene reference.

Saving the Debye plot

The parameters and data inputted to produce the plot can be saved by pressing the **Save** button, and then reviewed at a later stage by pressing the **Load** button.

Copying the Debye plot

The graph can be pasted into another application (such as Microsoft Word or Excel) by selecting the **Copy** button.

Results area

With both the table data and sample parameters entered the results will be automatically calculated and shown alongside the graph. The results displayed are:

- **MW (kDa) - molecular weight**

Shows the measured weight of a molecule within the sample expressed in atomic mass units; indicated in kiloDaltons. It is calculated from the intercept of the KC/RoP vs concentration Debye plot.

- **+/-**

The expected error in the molecular weight, derived from the scatter in the data about the least squares best fit line.

- **% Error in MW**

The error in the calculated molecular weight arising from the use of only a single angle. For isotropic scattering particles (diameter < ~50nm), this error should be negligible.

- **A₂ (ml mol/g²) - second virial coefficient**

A property describing the interaction strength between the molecule and the solvent. This is calculated from the slope of the plot.

- **K (Mol cm²/g)**

The instrument optical constant.

Molecular weight (MW) and shape estimates

The hydrodynamic size measured by Dynamic Light Scattering (DLS) is defined as “the size of a hypothetical hard sphere that diffuses in the same fashion as that of the particle being measured”. In practice though, macromolecules in solution are non-spherical, dynamic (tumbling), and solvated. Because of this, the diameter calculated from the diffusional properties of the particle will be indicative of the apparent size of the dynamic hydrated/solvated particle. Hence the terminology, **‘Hydrodynamic’ diameter**.

If the **molecular weight** (or mass) and the partial **specific volume** (inverse density) for the particle being measured are known, then a mass equivalent spherical size can be calculated. The closer the particle is to being spherical, the closer the mass equivalent spherical diameter will be to the DLS measured hydrodynamic diameter. In fact, it is the difference in these two values, coupled with **Perrin** theory (below), that allows **particle shape** information to be extracted from DLS measurements.

Once the **molecular weight** - either measured or estimated - and the **specific volume** are known, the **particle shape** information can be estimated by using the **Protein utilities** tool. The Shape estimate calculator takes the entered data and then applies two equations - the **Stokes-Einstein** and the **Perrin** factor.

Stokes-Einstein equation

The measured data in a DLS experiment is the correlation curve. Embodied within this curve is all of the information regarding the diffusion of particles within the sample that has been measured.

By fitting the correlation curve to an exponential function, the diffusion coefficient (D) can be calculated (D is proportional to the lifetime of the exponential decay).

With the diffusion coefficient (D) now known, the **Hydrodynamic diameter** can be calculated by using a variation of the Stokes-Einstein equation.

The **Stokes-Einstein** equation for the Hydrodynamic diameter is :

$$D_H = \frac{kT}{f} = \frac{kT}{3\pi\eta D}$$

D_H : Hydrodynamic diameter.

k : Boltzmann constant.

f : Particle frictional coefficient.

η : Solvent viscosity.

T : Absolute temperature.

D : Diffusion coefficient.

The Stokes-Einstein equation was developed using the assumption of hypothetical hard spheres.

Perrin factor

For non-spherical particles, the **Perrin** or shape factor (F) can be used to estimate particle shape.

The Perrin factor is used to calculate the prolate and oblate axial ratios for ellipsoids with the same Perrin factor value.

The Perrin factor is defined as the ratio of the frictional coefficient for a sphere with the same volume as the particle being measured, to the frictional coefficient for a sphere with the same mass as the particle being measured.

The **Perrin** factor (F) is :

$$F = \frac{f_{Vol}}{f_{Mass}} = \frac{6\pi\eta D_{Vol}}{6\pi\eta D_{Mass}} = \frac{D_{Vol}}{D_{Mass}} = \frac{D_H}{D_{Mass}}$$

D_H : Hydrodynamic diameter. The diameter as measured via DLS.

D_{Mass} : The diameter by mass. This is calculated from the known molecular weight and the specific volume of the particle.

f : Particle frictional coefficient.

η : Solvent viscosity.

T : Absolute temperature.

Shape estimate calculation

Enter the **Molecular weight** result, **Specific volume** and **Hydrodynamic diameter** (measured using dynamic light scattering) into the appropriate text boxes.

The Perrin (shape) factor (F), plus the Prolate and Oblate axial ratio will be automatically calculated and displayed in the results area.

If required, a solvent layer can be subtracted from the hydrodynamic radius when calculating the Perrin factor. To do this select the **Subtract solvent layer** check box.

Molecular weight estimate


Whilst the preferred method of measuring absolute molecular weight is by performing concentration dependent light scattering measurements, this can be very time consuming, from the point of view of the sample preparation. Providing that only an estimate of the molecular weight is required then it can be derived by utilising a relationship between the Hydrodynamic diameter and the molecule conformation.

To find the molecular weight estimate, enter the measured **Hydrodynamic diameter** value into the text box and the estimated molecular weight will automatically be calculated. The molecular weight is displayed in four ways:


- Globular
Proteins**
- Linear
polysaccharides**
- Branched poly-
saccharides**
- Starburst
polymers**

Hydrodynamic diameter estimate

This works in the same way as above, except the **molecular weight** (in kDaltons) is entered to find the Hydrodynamic diameter instead.



Concentration utilities



Select the **Concentration utilities** tab to view concentration and scattering parameters.

Concentration and Scattering

This area of the window contains features to establish the concentration and scattering levels that may be observed from the sample.

Enter the values from the measurement into the table. On entering each value, press the return key afterwards and the results table will be updated.


Minimum Concentration Calculator

By entering only the molecular weight the Sample concentration (mg/ml) values required for performing a measurement can be calculated:

- Sizing sample concentration on a Zetasizer Nano **S** instrument.
- Sizing sample concentration on a Zetasizer Nano **S90** instrument.
- Zeta potential sample concentration on a Zetasizer Nano **Z** instrument.

Scattering functions

A scattering function plot can be generated by inputting the measurement data in the list on the right, in this window:




Enter the values from the measurement into the table and press **Plot** - the graph will update to show the values entered.

To see the result of changing a value, change the required value and press **Add plot** - a new plot line will be added to the graph.

Protein charge & $f(Ka)$

This tool is used to calculate $f(Ka)$ for the Henry equation using known size and ionic strength. This tool uses the Ohshima equation for monovalent salts. The second calculation calculates protein charge from the electrophoretic mobility and Stokes radius.



Henry equation

The basic zeta potential measurement performed by the Zetasizer Nano measured the electrophoretic mobility of the particle or molecule under investigation. The mobility is related to the zeta potential using the Henry equation where:

$$U_E = \frac{2\varepsilon z f(\kappa a)}{3\eta}$$

where:

U_E = electrophoretic mobility

Z = zeta potential

ϵ = dielectric constant

η = viscosity

$F(\kappa a)$ = Henry's function.

- Henry's function ($f(\kappa a)$) is defined as :

$$\frac{1}{k} (nm) = [1.989 \times 10^{-3}] \left[\frac{\epsilon_r T(K)}{I(M)} \right]^{1/2}$$

where:

ϵ_r = dielectric constant

ϵ_0 = permittivity of free space (8.8542×10^{-12} C/Vm)

a = hydrodynamic radius

κ = inverse Debye length (1/nm)

NA = Avogadro's number (6.022×10^{23} /mole)

K = Boltzmann's constant (1.38065×10^{-23} m²kg/s²K)

T = temperature (K)

I = ionic strength (moles/L)

This calculator allows a more specific value other than the Smoluchowski or Huckel estimates to be calculated and used for zeta potential measurements.

Protein charge

The protein charge calculation calculates protein charge from the measured electrophoretic mobility and the hydrodynamic size. The charge is calculated from the following equation:

$$Z = [5.986 \times 10^{-5}] \epsilon_r T(K) \kappa a^2 (nm) \times \left[2 \sinh \left(\frac{5801.4 \zeta(V)}{T(K)} \right) + \frac{4}{\kappa a} \tanh \left(\frac{2900.7 \zeta(V)}{T(K)} \right) \right]$$

where:


Z = calculated protein charge

ζ = the zeta potential

and the other values can be taken from the above.

Virial Diameter

A tool to calculate the 'virial diameter' from the measured molecular weight and 2nd virial coefficient (A₂).




Dynamic (DLS) Debye plot

A dynamic Debye plot can be created by inputting sample parameters and data rather than using measured data. This feature can be useful for various reasons, for example:

- A dynamic Debye plot can be created by inputting the concentration points from an existing measurement, then additional concentration points can be added.
- Any of the parameters in the dialogue box can be changed, and the other parameters are instantly recalculated. This can be used to investigate the sensitivity of the result to changes in any parameter. For example, by first entering the concentration points from an existing measurement, a sample parameter - e.g. sample temperature - can be altered, with the effect immediately observed on the Debye plot. This saves time in performing the original measurement again at the different temperature.

The format of the plot can be altered by choosing **Chart properties** from the drop-down menu in the top right hand corner of the graph.



Adding and editing sample parameters and table data

To access the **Debye plot** select **Tools-Calculators** and then the **Debye plot** tab. Select the **Table** tab to begin entering data into the table.

- To **define a new** concentration, press the **Add...** button. The table **Input values** window will be displayed
 - Specify the **Concentration** and **Measured diffusion coefficient** values - either new values or ones taken from an existing measurement.
- To **modify** a concentration, select it from the list and press the **Modify...** button. The **Input values** window is displayed, allowing the parameters to be changed.
- A concentration can be **deleted** by selecting the concentration from the list and pressing the **Delete...** button.
- Select the **Graph** tab to see the resultant **Debye plot**.
- The table values and graph plot can be subsequently altered by changing the **Sample, Data** and **System** parameters in the measurement parameters table on the right of the window. These parameters are described in the following section.

Measurement parameters table

When all the concentration values have been added into the table, use the measurement parameters table on the right hand side to alter the result and dynamic Debye plot.

The parameters are:

- **Molecular weight**
The molecular weight of the sample (optional).
- **Partial specific volume**
The partial specific volume of the sample (optional).
- **Second virial coefficient**
The second virial coefficient of the sample (optional).
- **Solution Viscosity**
The viscosity of the solution.
- **Temperature**
The measurement temperature. This is displayed in both Celsius and Kelvin. If either is completed, the other will be updated.

Saving the dynamic Debye plot

Click **Save** to save the parameters and other inputted data. To reviewed the saved parameters at a later stage, click the **Load** button.

Copying the dynamic Debye plot

To paste the graph into another application (such as Microsoft Word or Excel) select the **Copy** button.

Results area

When both the table data and sample parameters have been entered the results are automatically calculated and shown alongside the graph. The results displayed are:

- **Diffusion coefficient**
The calculated diffusion coefficient at zero concentration.
- **Stokes radius**
The calculated radius at zero concentration.
- **DLS interaction parameter**
The DLS interaction parameter (kD) is also known as the dynamic virial coefficient and is related to the slope of the line on this plot.
- **Frictional coefficient**
The frictional coefficient.

Interparticle distance

This tool is used to calculate the distance between the particles based on their concentration and molecular weight. It also estimates the thickness of the electrostatic layer based on protein charge and ionic strength.

The values are calculated according to the following formulae:

- **Dielectric constant** of the media ϵ_r is:

$$\epsilon_r = (-1.204 \times 10^{-6})T^3 + (1.879 \times 10^{-3})T^2 - (1.162)T + (289.9)$$

- **Hydrodynamic radius** estimated from molecular weight based on a globular protein model:

$$Rh = 0.294 \times Mw^{0.428}$$

- **The Debye length** is:

$$\frac{1}{k} (nm) = [1.989 \times 10^{-3}] \left[\frac{\epsilon_r T(K)}{I(M)} \right]^{1/2}$$

- **Separation distance** is:

$$Separation\ distance = (0.001C)^{1/3} \times Rh$$

where:

ϵ_0 = permittivity of free space (8.8542×10^{-12} C/Vm)

κ = inverse Debye length (1/nm)

N_A = Avogadro's number (6.022×10^{23} /mole)

K = Boltzmann's constant (1.38065×10^{-23} m²kg/s²K)

T = temperature (K)

I = ionic strength (moles/L)

This assumes that the available solvent volume is equally distributed between all particles, and that the particle is located in the centre of a cube that would represent this solvent volume 'occupied' per particle. The interparticle distance, on average, is then given by the edge length of the cube, representing the distance from the centre of one cube to the centre of the neighbouring cube.

Mixture Viscosity

A simple tool that calculates the overall viscosity of a mixture of solutions based on a volume weighted mean viscosity.

Index

A

Acceptable Zeta Difference 5-6
Advanced features
 Flow-mode 4-1

B

Baseline pointer 4-7
Branched polysaccharides 6-19

C

Calculators tool 6-1
Cell
 Choosing 2-3
 Cleaning 2-7, 2-13
 Filling 2-7, 2-13
 Triangle 2-22
Chromatogram summary 4-8
Chromatographic column 4-1
Concentration 6-20
Concentration utilities 6-10, 6-19
Cuvette holder 2-2, 2-23
Cuvettes 2-3

D

Debye plot 6-12
Delay volume 4-11
Disposable cuvettes 2-3, 2-4
Disposable low volume polystyrene cuvette 2-5
dn/dc 6-15

E

Enable External input 4-10
Exporting the Flow-mode 4-4, 5-15
External Input 1 and 2 (SOP) 4-10

F

Filling the cell 2-7, 2-13
Flow settings (SOP) 4-13
Flowcell 2-10
Flow-mode 4-1, 4-9

 Connectivity 4-3
 Exporting the data 4-4
 Process monitor 4-3
Flow-mode (Editing) 4-5
Flow-mode (Inspecting) 4-5
Flow-mode measurements 4-14
Flow-mode measurements - Standard report 4-15, 6-9
Flow-mode option 1-4
Folded capillary cell 1-2, 2-3, 2-4
 Insertion 2-16
 Orientation 2-17

G

Gain 4-11
Glass - round aperture cuvette 2-5
Glass - square aperture cuvette 2-5
Glass cuvettes 2-3
Globular Proteins 6-19

H

High Concentration cell 2-25, 2-29
High concentration cell 2-25
High temperature 1-4
HT label 1-4
Hydrodynamic diameter 6-10, 6-16
 Estimate 6-19
Hydrodynamic radius 6-15

L

Linear polysaccharides 6-19
Low volume cell
 Filling 2-9
Low volume Glass flow cuvette 2-6
Low volume quartz cuvette 2-6

M

Measurement (SOP) 4-12
Measurement tab 4-8
Microrheology 5-1
 SOP 5-4

Minimum Concentration Calculator 6-20
Model number 2-2
Molecular weight
 Calculations 6-16
 Estimate 6-18
Molecular weight calculations 6-10
MPT-2 Autotitrator 1-4

N

Narrow band filter 1-5

O

Offset 4-11
Option label 2-2

P

Peaks 4-6
Peaks tab 4-8
Perform size measurement 6-6
Perrin factor 6-17
Process monitor 4-3

Q

Quartz cuvettes 2-3

R

Refractive index 6-15
Report
 Flow-mode measurements 4-15, 6-9
Reports 4-15, 5-13, 6-9
Result tab
 Flow-mode 4-14

S

Scattering 6-20
Scattering functions 6-10
Second virial coefficient 6-16

Serial number 2-2
Shape 6-10
Shape estimate calculation 6-18
Shape model 6-15
Size measurement (SOP) 6-6
SLS Debye plot 6-10
SOP dialogues
 External Input 1 and 2 (Flow-mode) 4-10
 Flow settings (Flow-mode) 4-13
 Measurement (Flow-mode) 4-12
Starburst polymers 6-19
Stokes-Einstein equation 6-17

T

Thermal cap 2-9
Thermal contact plates 2-2, 2-16
Tools 4-6
Tracer compatibility 5-6
Tracer concentration 5-7
Tracer settings 5-5
Traces 4-6
Traces tab 4-8

U

Universal 'Dip' cell 2-3
 Filling 2-20
 Insertion 2-21
 Removal 2-22

V

Viewing tools 4-6, 4-7

W

Weld line 2-18

Z

Z-average mean 4-5